Reference Device and Lok-Bar™

REFERENCE GUIDE

CIVCO

CE

Towards Gantry

Towards Gantry
Směrem k nosné konstrukci
Mod ophængning
Naar rijbrug
Portaali suunas
Pöytäa kohti
Couvercle de portique
Zum Gestell hin
Προς πλαίσιο
A gantry irányába
Verso il gantry
構台方向
Mot gantry
W kierunku suwnicy
Em direção à estrutura
Em direção à torre
朝向扫描架
Spre suportul de scanare
К гентри
Smerom k portálu
Hacia la pasarela
Mot gantry
Kızağa Doğru

INTENDED USE

The device is intended to support and aid in positioning a patient during radiologic and other medical procedures.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.



⚠ WARNING

- Do not use if device appears damaged.
- Regularly inspect o-rings to ensure device fits properly.
- Do not use to secure positioning devices.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Only use CIVCO compatible accessories.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.


USING REFERENCE LOK-BAR™

1. Snap Reference Lok-Bar™ at appropriate longitudinal position.
2. Use reference marks on Reference Lok-Bar™ to determine position of system in reference to couchtop.

USING REFERENCE DEVICE

1. Attach Lok-Bar™ to couchtop at 0 indexing location on side of couchtop.
2. Place reference device on Lok-Bar™ with arrows pointing towards gantry.
3. Use marks on reference device to determine position of couchtop relative to system.

REPROCESSING**WARNING**

Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Disinfect with a mild, non-abrasive cleanser.
2. Wipe dry with lint-free cloth or allow to air dry prior to reassembly and use.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

STORAGE

NOTE: Device should be stored to keep all components together and secure.

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COPYRIGHT © 2020 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. LOK-BAR AND PRODIGY ARE TRADEMARKS OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.