Headrests

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INTENDED USE

The device is intended to aid in patient positioning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

- *Verify all angles of treatment and attenuation characteristics prior to treating patients.*
- *Product shape may have slight variations; the same headrest must be used for a patient throughout setup and treatment cycles.*
- *Do not use if device appears damaged.*
- *Ensure accessory device is mounted accurately on headrest prior to patient treatment.*
- *Ensure device is secure prior to use.*
- *Verification image should be taken to confirm patient position prior to treatment.*
- *Do not reposition device with patient on it.*
- *Ensure patient remains still for duration of setup and treatment.*

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

USING HEADRESTS

- NOTE:
- **Timo:** headrest must have CQ Medical interface plate attached to underside to use with baseplate.
 - **Silverman:** headrest must have indexing holes on lateral edges to use with baseplate or Type-S™.
 - **Timo and Silverman** are not for use with Posifix® baseplates.
 - **Posifix® headrests:** If desired for use with Type-S™, Posi-S™ adapter is needed. If desired for use with Uni-frame® baseplate, cushion adapter is needed. Posi-S™ adapter is not MR environment safe.
1. Place headrest onto table top, into recessed area of Uni-frame® supine baseplate, or into recessed area on Posifix® baseplate.
 2. Place lowest elevation of the headrest at superior end of table allowing curved surface to properly position the neck.

NOTE: 301005BLK/301205 headrest must be positioned with highest elevation of headrest at superior end of table.

3. Position patient on headrest.

REPROCESSING

WARNING

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*

1. If necessary, remove gross visible contaminants with warm soapy water and a soft cloth.
2. Disinfect with a mild, non-abrasive cleanser.

MAINTENANCE

NOTE: Inspect for damage such as cracks, breakage, or warpage.

STORAGE

NOTE: Store headrests on a flat surface. Do not stack.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. TYPE-S IS A TRADEMARK OF MEDTEC LLC. POSIFIX AND UNI-FRAME ARE REGISTERED TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MEDETEC LLC