Vac-Lok™ Cushions

INSTRUCTIONS FOR USE (IFU)

CE

INTENDED USE

The device is indicated to aid in supporting and positioning adult and pediatric patients undergoing radiation therapy including electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

- Do not use if device appears damaged.
- Verify all angles of treatment, attenuation characteristics and WET values prior to treating patients. Refer to Technical Data Sheet at www.CQmedical.com for WET values.
- Do not reposition device with patient on it.
- Remove sharp objects from patient and treatment area to avoid puncturing.
- Ensure Vac-Lok™ bag is deflated to proper vacuum prior to each use.
- Never leave patient unattended on flaccid (inflated) Vac-Lok™ cushion.
- Do not use Lok-Bar™ to secure positioning devices.
- Ensure cushion is not over-inflated to avoid bursting or difficulty in positioning patient.
- Nylon fabric Vac-Lok™ cushions may become soiled if not covered with sheet.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

MRI SAFETY INFORMATION

MR
Conditional

- Non-clinical testing and scientific rationale has demonstrated that the Vac-Lok™ Cushion is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:
 - Static magnetic field of 1.5 T and 3.0 T
 - Maximum spatial field gradient of 2,000 gauss/cm (20 T/m)
 - Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Operating Mode)
- Under the scan conditions defined above, the device is expected to produce a maximum temperature rise of less than 0.4°C after 15 minutes of continuous scanning.
- In non-clinical testing, the image artifact caused by the device extends up to 82mm from the Vac-Lok™ valve with stem when imaged with a spin echo pulse sequence and 3.0 T MRI System.

NOTE:

- Vac-Lok™ cushion should be flaccid prior to positioning patient on it.
- Ensure integrity of Vac-Lok™ cushion is not compromised during daily handling.
- Place female valve connector onto cushion valve to let air into Vac-Lok™ cushion without using pump.

USING VAC-LOK™ CUSHIONS

1. Place appropriate Vac-Lok™ cushion onto tabletop.
 - Vac-Lok™ F-Series cushion: place onto Lok-Bar™ pins (*refer to Lok-Bar™ instruction*).
 - Cushions with multiple indexing points: attach cushion with one fixed locking feature and one Slide Guide Lok-Bar™.
2. Smooth Vac-Lok™ cushion to evenly disperse beads.
3. Attach vacuum hose to cushion valve.

WARNING

- Locate valve away from region of interest to prevent image artifacts.

4. Place sheet over cushion.
5. Position patient onto Vac-Lok™ cushion. Distribute beads around anatomical regions to be immobilized.
6. Evacuate air from Vac-Lok™ cushion while forming around anatomical regions until desired level of vacuum is achieved.

REPROCESSING

WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Remove visible contaminants using a soft damp cloth.
2. Wipe dry with lint-free cloth or allow to air dry.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. VAC-LOK AND LOK-BAR ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
**MEDTEC LLC**