# Hipstop

---

REFERENCE GUIDE

**CIVCO**

CE

**INTENDED USE**

Device is intended to aid in positioning patients on a breastboard.

**CAUTION**

Federal (United States) law restricts this device to sale by or on the order of a physician.


**⚠ WARNING**

- Do not use if device appears damaged.
- Ensure patient does not shift during treatment.
- Ensure device is attached securely prior to treatment.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at [WWW.CIVCORT.COM](http://WWW.CIVCORT.COM).
- Verify patient position with completed setup sheet prior to treatment.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Ensure all accessories are attached securely prior to treatment.
- Do not reposition device with patient on it.
- Do not allow patient to reposition themselves.


**NOTE:** If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

**ASSEMBLING HIPSTOPS****MT-250-100:**

1. Lift transition pad.
2. Remove six (6) screws and attachment bar from breastboard.


3. Center hipstop with breastboard and using existing screws insert one (1) onto each end of the new attachment bar.
4. Lift hipstop and insert remaining four (4) screws into attachment bar.

**MT-350-300:**


1. Remove existing hinge plate from breastboard.
2. Attach hipstop with hingeplate and screws (*provided*).

3. Ensure cap screw is attached to hipstop link.

MT-350-302:

1. Remove existing hipstop.
2. Align brackets on hipstop with breastboard.
3. Attach hipstop with existing screws (2).


MT-350-303:

**WARNING**

- Do not use Lok-Bar™ to secure positioning devices.

1. Remove existing hipstop from breastboard.
2. Place hipstop onto Lok-Bar™ pins (Refer to Lok-Bar™ instruction).

REPROCESSING**WARNING**

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*

1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2020 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. LOK-BAR IS A TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.