
Iso-Align™

INSTRUCTIONS FOR USE (IFU)

I. Side Mounts	VIII. Faceplate
II. Leveling Feet (2)	IX. Retaining Clamp
III. Baseplate	X. Knurled Screw
IV. Baseplate Crosshairs	XI. Iso-Point
V. ODI (Optical Distance Indicator) Checking Tool	XII. Lok-Bar™ Two Pin Mount
VI. Tungsten Rod	XIII. Faceplate Circle
VII. Film Slot	XIV. Tolerance Indicators

INTENDED USE

The device is intended to provide secondary radiologic quality assurance of linear accelerators and room lasers following completion of primary system calibration.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING


- Do not use if device appears damaged.
- Refer to your system manual for device alignment if any test fail.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

ISO-ALIGN™ SETUP

BASELINE CHECK

1. Ensure collimator is square with gantry rotation.
 2. Secure tabletop and position at 0° (180° Varian convention).
- NOTE: Ensure gantry and collimator rotations do not interfere with Iso-Align™ during planned operations.
3. Place Iso-Align™ onto Lok-Bar™ pins (refer to Lok-Bar™ instruction) or directly onto tabletop.
 4. Align collimator and baseplate cross hairs.
 5. Level baseplate with digital level (provided) by adjusting leveling feet.
 6. Insert faceplate into baseplate.
 7. Secure faceplate and Iso-Point to baseplate with retaining clamps and knurl screws at side mounts.
- NOTE: Retaining clamp is not used on side with Iso-Point. Use Knurled screw only.

8. Check collimator and baseplate cross hairs alignment.

VERIFY ISOCENTER

1. Firmly attach a mechanical pointer to collimator.

WARNING

- Ensure Iso-Point tip and collimator pointer do not come in contact at any point in rotation.

NOTE: Collimator should read 100cm for a 100cm SAD (source to axis distance) accelerator.

2. Adjust tabletop longitudinally and vertically until pointer and Iso-Point tips are closely aligned.
3. Rotate gantry. Note position of collimator pointer compared to Iso-Point at 0°, 90°, 180°, and 270° and additional desired angles.

NOTE: If collimator pointer tip stays within Iso-Point circle, mechanical isocenter of gantry is +/- 1.5mm or less. Repeat process if desired results are not achieved.

4. Bring gantry to vertical.
5. With faceplate horizontal, move table toward gantry until collimator and faceplate cross hairs align.
6. Check distance to faceplate with collimator pointer.
7. Ensure faceplate is at isocenter. Rotate faceplate to vertical.
8. Rotate gantry clockwise 90°. Check collimator and faceplate cross hairs alignment.
9. Return gantry to vertical.

USING ISO-ALIGN™

LASER ALIGNMENT

1. Rotate faceplate to vertical.
2. Check vertical and horizontal alignment of lateral and ceiling lasers.
3. Check back pointer laser alignment with end plate scribe line.
4. Check sagittal laser alignment with opposite end plate scribe line.

SYMMETRIC LIGHT FIELD ACCURACY

1. Rotate faceplate to horizontal.
2. Open upper jaws to maximize field size.
3. Close lower jaws so light field edge is aligned with 5cm scribe line. Check for variance.
4. Repeat for 10cm, 15cm and 20cm.

NOTE: Use tolerance indicators located at corners of faceplate fields to verify light field falls within +/- 1mm of scribe line.

ASYMMETRIC LIGHT FIELD ACCURACY


1. Switch to asymmetric mode (if applicable).
2. Set one asymmetric collimator to 10cm. Check opposite jaw.
3. Move opposite asymmetric jaw to 2.5cm. Record digital readout. Verify and record reading for jaw and total field size.
4. Repeat for 5cm, 7.5cm and 10cm.
5. Open lower jaws to maximize field size. Repeat steps 2–4 for upper jaws.
6. Calculate differences between actual field size and readouts. Check for agreement with your standard tolerances.

NOTE: Use tolerance indicators located at corners of faceplate fields to verify light field falls within +/- 1mm of scribe line.

LIGHT AND RADIATION FIELD COINCIDENCE

1. Rotate faceplate to horizontal.
2. Place 25.4cm x 30.5cm Ready-pack™ film into faceplate film slot.
3. Set appropriate field size.
4. Expose film.

NOTE: • Light field should align with tungsten rods on Iso-Align™.
• For 15 x 15cm field, tungsten rods are located at 1.6mm and 3.2mm +/- 0.2mm from scribe line as shown below.

5. Repeat for other appropriate field sizes and energies.

COLLIMATOR MECHANICAL ISOCENTER, CROSS HAIRS ISOCENTER, AND COLLIMATOR READOUT ACCURACY

1. Place mechanical pointer in appropriate accessory mount and position pointer tip at faceplate center.
2. Rotate collimator. Pointer tip should remain within faceplate circle.

NOTE: Diameter of faceplate circle is 3mm +/- 0.2mm.

3. Remove pointer and rotate collimator. Cross hairs center should remain within faceplate circle.
4. Use corners of faceplate fields for accurate 45° settings and faceplate crosshairs for 90° settings.

COLLIMATOR RADIATION ISOCENTER


1. Place 25.4cm x 30.5cm Ready-pack™ film into faceplate film slot.
2. Open both jaws to 25cm or appropriate field size. Expose film.
3. Open lower jaws fully.
4. Close upper jaws to approximately 1mm. Expose film.
5. Rotate collimator preferred increments. Expose film.
6. Check for pattern variance.
7. Repeat procedure for lower jaws.

CEILING LASER ALIGNMENT, GANTRY ANGLE VERIFICATION, CROSS HAIR LOCATION WITH GANTRY ANGLE, LIGHT VS. RADIATION FIELD COINCIDENCE AT OTHER GANTRY ANGLES

1. Rotate gantry. Observe ceiling laser position on horizontal faceplate. Observe laser back pointer location on end plate as gantry is rotated.
2. Position faceplate at appropriate angle.
 - NOTE:
 - Detent positions are provided every 10 degrees ± 0.4 degrees (0.175 radians ± 0.007 radians).
 - Cross hairs center should remain within faceplate circle.
 - Diameter of faceplate circle is $3\text{mm} \pm 0.2\text{mm}$.
3. Compare light field with faceplate markings.
4. Verify gantry angle.
5. Verify gantry digital readouts and readouts on record. Verify system.
6. Verify collimator pointer readings at 90° and 270° .
7. Take radiation field vs. light field films.

OPTICAL DISTANCE INDICATOR (ODI) ACCURACY

1. Place gantry at 0° .

NOTE: ODI block tolerance is $\pm 0.6\text{mm}$.

TREATMENT TABLE ISOCENTER VERTICAL AND HORIZONTAL TRAVEL


1. Rotate gantry to vertical.
2. Rotate faceplate to horizontal.
3. Check collimator and faceplate cross hairs alignment.
4. Rotate table. Collimator cross hairs should remain within faceplate circle.
 - NOTE: Diameter of faceplate circle is $3\text{mm} \pm 0.2\text{mm}$.
5. Verify table rotation angle reading.
 - NOTE: Use faceplate field corners for accurate 45° settings and faceplate cross hairs for 90° settings.
6. Return couch to 0° .
7. Check collimator and faceplate cross hairs alignment.
8. Move couch vertically. Collimator cross hairs should remain within faceplate circle.
9. Rotate faceplate to vertical.
10. Move table in and out horizontally. Horizontal side laser cross hairs should remain stable with faceplate cross hairs.

RADIATION ISOCENTER

1. Rotate Iso-Align™ to horizontal.
2. Level device with digital level.
3. Align side lasers with cross hairs on each end plate.
4. Verify faceplate top is at isocenter.
5. Check collimator and faceplate cross hairs alignment.
6. Place $25.4\text{cm} \times 30.5\text{cm}$ Ready-pack™ film into faceplate film slot.
7. Open both jaws to 25cm or appropriate field size. Expose film.
8. Rotate faceplate to vertical. Verify sagittal laser is at a true vertical position.
9. Close collimator jaws perpendicular to long axis of film to approximately 1mm .
10. Close other collimator to 5cm . Expose film.
11. Repeat at appropriate gantry angles.

STORAGE

- NOTE:
- Place end cap (*provided*) over Iso-Point tip before storage.
 - Store device in original packaging.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. ISO-ALIGN AND LOK-BAR ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS.
**MEDTEC LLC**