Waterbath

REFERENCE GUIDE

CIVCO
Radiotherapy™

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Certified to IEC 61010-1:2010 Ed. 3+A1
Certified to IEC 61010-2-010:2019 Ed. 4
Certified to IEC 60601-1-2:2014 Ed. 4
Conformed to CSA C22.2 61010-2-010:2019 Ed. 4

Device Specifications	REF# MTPVEXA, MTPVEXC, MTPVEXD, MTPVEXG, MTPVEXH	REF# MTPV1
Classification	Class 1	Class 1
Operating Voltage	120 VAC	230 VAC
Current (Max Load)	10 Amps	5 Amps
Frequency	60 Hz	50 Hz
Fuse	15 Amp (Slow Blow)	10 Amps (Fast Blow)
Weight	~ 70 lbs	~ 70 lbs

Replacement Parts - Contact CIVCO Radiotherapy for service.	
MTRH1500	Power Cord EJCFCSETK3
MTRH1500EXA	Power Cord EJCFCSETA
MTRH1500EXC	Power Cord EJCFCSETC
MTRH1500EXD	Power Cord EJCFCSETD
MTRH1500EXG	Power Cord EJCFCSETG
MTRH1500EXH	Power Cord EJCFCSETH
MTRH1500EXI	Power Cord EJCFCSETI
MTRH1500RK	120V Inlet Repair Kit: Inlet, Switch, Power Cord, Wires and Terminals
MTRHF15	Fuse
MTCGPLLG	Pan liner
MTRHGS01	Gas Strut

Test / Standard	Emissions Class and Group / Immunity Test Level
CISPR 11, Conducted Emissions	Class B; Group 1
CISPR 11, Radiated Emissions	Class B; Group 1
IEC 61000-3-2, Harmonic Current Emissions	Class A
IEC 61000-3-3, Voltage Changes, Voltage Fluctuations, and Flicker	All Parameters
IEC 61000-4-2, Electrostatic Discharge Immunity	±8kV contact ±2kV, ±4kV, ±8kV, ±15kV air
IEC 61000-4-3, Radiated RF Immunity	3V/m, 80-2700MHz, 80% 1kHz AM
IEC 61000-4-3, Proximity fields from RF wireless equipment	Section 8.10 of the IEC 60601-1-2 standard
IEC 61000-4-4, Electrical Fast Transient / Burst Immunity	2kV, 100kHz repetition rate
IEC 61000-4-5, Surge Immunity	±0.5kV, ±1kV Line-to-Line ±0.5kV, ±1kV, ±2kV Line-to-Ground
IEC 61000-4-6, Immunity to Conducted Disturbances, Induced by Radio Frequency Fields	3V, 0.15-80MHz, 80% 1kHz AM 6V in ISM Band within 0.15-80MHz, 80% 1kHz AM
IEC 61000-4-8, Magnetic Field Immunity	30A/m, 50 or 60Hz
IEC 61000-4-11, Voltage Dips	0% (100% reduction), 0.5 cycle 0% (100% reduction), 1 cycle 70% (30% reduction), 0.5 sec
IEC 61000-4-11, Voltage Interruptions	0% (100% reduction), 5 sec

PRIOR TO USE:

- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally. Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the waterbath, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Waterbath heats and maintains temperature of water, allowing thermoplastic material to be heated. If performance is lost, due to electromagnetic disturbances, thermoplastic mask may not be properly heated.
- Failure to use the specified fuse could result in device damage or injury.
- The use of accessory equipment and/or hardware not complying with the equivalent product safety and EMC requirements of this product may lead to a reduced level of safety and/or EMC performance of the resulting system. Consideration relating to the choice of accessory equipment used with this product shall include:
 - The use of the accessory in the patient vicinity.
 - Evidence the safety certification of the accessory has been performed in accordance with the appropriate IEC 60601-1 Harmonized National Standards.
 - Evidence the EMC certification of the accessory has been performed in accordance to the IEC 60601-1-2 Harmonized National Standards.
- Two person lift required when lifting waterbath. Single person lift could cause injury. Use assistance and proper lifting techniques when moving or replacing. Waterbath weight is 32 kg (70 lbs) when empty. Lift only when waterbath is empty and cool. Lift by placing hands under waterbath. Do not use handle to life entire waterbath or harm could occur.
- Prior to use, equipment must be assembled by inserting one end of the power cord in the waterbath and connecting the other plug to a GFCI circuit rated for your waterbath. Use only power cords provided by CIVCO Radiotherapy. Failure to comply could result in harm.
- Equipment shall be operated indoors only and at temperatures between 4-41°C (40-105°F).
- Equipment shall operated at altitudes of 2000 meters (6561 ft) or less.
- Equipment shall be operated at relative humidity of 70% or less.
- Equipment shall be operated within 2% of the rated input voltage.
- Only non-conductive pollution occurs except occasionally a temporary conductivity caused by condensation is normal.
- Equipment shall be mounted only on surfaces capable of safely bearing at least 284 kg (625 lbs).
- Equipment shall be operated with CIVCO Radiotherapy supplied power cord plugs that match the wall socket. Use of adapters or transformers is not allowed.
- Equipment shall be operated with medical grade ventilation.
- All sides of waterbath shall be positioned at least 10 cm (4 inches) from a wall.
- Waterbath liners can be changed often, as desired.
- Waterbaths shall be returned to CIVCO Radiotherapy for servicing by authorized personnel. Attempting to self-service could result in harm.

INTENDED USE

Device is intended to heat water.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Use only power cords provided.
- Do not use if device appears damaged.
- Hot water poses a scalding hazard.
- To prevent pan liner from melting, periodically check waterbath to ensure water level does not drop below 25mm (1") from bottom of waterbath.
- Routinely check waterbath and power cords for signs of deterioration. Do not use if device appears damaged.
- Verify all settings and temperatures prior to using device.
- If equipment is not used as specified, the protection provided by the equipment may be impaired.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

WATERBATH SETUP**⚠ WARNING**

- Ensure fan has sufficient ventilation and is not blocked.

1. Apply pipe thread compound to ball valve.
2. Attach ball valve to waterbath.
3. Plug waterbath into appropriate GFCI (*ground fault circuit interrupter*) receptacle.

⚠ WARNING

- Injury can occur if device is not properly grounded.

4. Remove all contents from waterbath tank and close lid.
5. Turn waterbath power on.
6. Set waterbath controls:
 - Tank Temperature: Press + or - as necessary to set tank temperature to approximately 75°C (165°F).
 - Selecting °C/°F unit: Simultaneously press + and - . Press + until 10 is displayed. Press + once more to toggle between Fahrenheit or Celsius option.
7. Allow waterbath to run at temperature settings for a minimum of 24 to 48 hours for thermal stabilization of heating element.

USING WATERBATH**⚠ WARNING**

- Touch only control panel and lid handle while waterbath is hot.

NOTE: Waterbath is suitable for use in Radiation Oncology planning and treatment environments.

1. Plug waterbath into appropriate GFCI (*ground fault circuit interrupter*) receptacle.
2. Place bath liner into waterbath to cover entire bottom surface.
3. Ensure waterbath contains 8cm (3 inches) of water. Close lid.

NOTE: • Do not use Reverse Osmosis (RO) water due to its corrosive nature. For optimal performance and to limit the buildup of deposits, use distilled water.
• Device is made of stainless steel. Do not put items containing ferric metal in tank.

4. Turn waterbath power on.
5. Set waterbath controls:
 - Tank Temperature: Press + or - as necessary to set tank temperature to approximately 75°C (165°F).
 - Selecting °C/°F unit: Simultaneously press + and - . Press + until 10 is displayed. Press + once more to toggle between Fahrenheit or Celsius option.
6. Allow waterbath two to three hours to heat water.

REPLACING FUSE


1. Turn off waterbath.
2. Unplug waterbath from electrical source.
3. Remove the cover from fuse holder.
4. Remove the fuse from fuse holder using a small flat head screwdriver.
5. Obtain correct fuse as listed in replacement parts table.
6. Insert fuse in fuse holder, securely seating fuse against terminals.
7. Place fuse holder cover back on the fuse holder.
8. Plug waterbath into electrical source.
9. Turn on waterbath. Power should be restored to waterbath.

REPLACING GAS STRUT

NOTE: • Ensure waterbath is turned off, unplugged, cooled and drained before performing maintenance.
• Ensure large end of gas strut is connected to waterbath lid.


MTRHGS01 / MTRHGS60

(Pronged end cap to fully stainless steel end cap)


MTRHGS01

(Polymer end cap with metal retaining clip to fully stainless steel end cap)

**MTRHGS01**

(Fully stainless steel end cap to fully stainless steel end cap)

**MAINTENANCE**

NOTE: Inspect device prior to use for signs of damage and general wear.

REPROCESSING**⚠ WARNING**

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

- NOTE:
- Drain and clean waterbath frequently. Failure to do so can cause contamination and deposit buildup.
 - Clean waterbath approximately once a week using a standard cleaning process for contamination control and removal of light buildup of deposits. More often cleaning may be needed depending on usage.
 - Clean waterbath approximately once a month with high grade stainless steel cleaner formulated to remove ferric metal from surfaces. More often cleaning may be needed depending on usage.
 - If unsure about compatibility of cleaning agent, contact CIVCO Radiotherapy.
- Turn waterbath off.
 - Allow water to cool to room temperature.
 - Attach hose to ball valve nozzle.
 - Turn ball valve nozzle to drain water from waterbath.
 - Clean inside of waterbath with appropriate cleaner. Rinse thoroughly.
 - Shut ball valve. Remove hose from ball valve nozzle.
 - Refill waterbath with new distilled water after each routine cleaning to help avoid cross contamination, ferric metal buildup and calcium buildup in tank.

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COPYRIGHT © 2022. ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. CIVCO RADIOTHERAPY AND CIVCO RT ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.