Proton Type-S™

For use with Hitachi Proton System

REFERENCE GUIDE

CIVCO
Radiotherapy™

CE

INTENDED USE

The device is indicated to aid in supporting, positioning, and/or immobilization of adult and pediatric patients undergoing radiation therapy of the head, brain, neck, and spine including radiosurgery and electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Do not use if device appears damaged.
- Verify all angles of treatment, attenuation characteristics and WET values prior to treating patients. Refer to Technical Data Sheet at WWW.CIVCORT.COM for WET values.
- Do not reposition device with patient on it.
- Ensure device is secure prior to use.


NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

USING PROTON TYPE-S™ OVERLAY

⚠ WARNING

- Do not use Lok-Bar™ to secure positioning devices.
- Do not extend inscribed line past end of table.

1. Place proton Type-S™ overlay onto Lok-Bar™ pins (refer to Lok-Bar™ instruction).
2. Wrap velcro straps around proton Type-S™ overlay and couchtop.


USING PROTON TYPE-S™ EXTENSION

⚠ WARNING

- Do not apply a load exceeding 308 lbs (140 kg).

1. Attach Proton Type-S™ Extension onto superior end of table.
2. Ensure Proton Type-S™ Extension is secured at pin and clamp attachments (refer to Hitachi System's instruction).

REPROCESSING

WARNING

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*

1. If necessary, clean surface by removing visual contamination with common germicidal wipe, such as alcohol. If visual contamination cannot be removed repeat cleaning steps and if necessary, discard device.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2022 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. CIVCO RADIOTHERAPY AND CIVCO RT ARE TRADEMARKS OF MEDTEC LLC. LOK-BAR AND TYPE-S ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.