# Alloy Melter

---

REFERENCE GUIDE

**CIVCO**

CE

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CAN/CSA C22.2 #1010.1-92 belgelidir

**INTENDED USE**

Device is intended to heat shielding alloy.

**CAUTION**

Federal (United States) law restricts this device to sale by or on the order of a physician.

**ALLOY MELTER SET-UP****⚠ WARNING**

- Do not use if device appears damaged.
- Consult your facilities manager before assembly and installation.
- Use only power cords provided.
- Do not move or slide alloy melter with molten alloy in tank.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

1. Close ball valve dispenser.
2. Plug alloy melter into appropriate GFCI (ground fault circuit interrupter) receptacle.

**⚠ WARNING**

- Ensure device is properly grounded to prevent injury.

3. Turn alloy melter on.

**ALLOY MELTER CONTROLS MUST BE SET IN THE FOLLOWING ORDER:**

4. Tank Temperature: Press + or - as necessary to set tank temperature to approximately 75°C (165°F).
5. Selecting °C/°F unit: Simultaneously press + and - . Press + until 10 is displayed. Press + once more to toggle between Fahrenheit or Celsius option.

NOTE: The procedure for changing the temperature units also adjusts the valve temperature set point. Step 6 must be performed after changing the temperature units to correctly set the valve temperature.

6. Set Valve Temperature: Simultaneously press + and - . Press + or - to set valve temperature setting. 3 for low temperature alloys and 5 for medium temperature alloys.
7. Remove alloy melter lid.
8. Allow alloy melter to run at temperature settings for a minimum of 24 to 48 hours for thermal stabilization of heating elements.

**USING ALLOY MELTER****⚠ WARNING**

- Routinely check alloy melter and power cords for signs of deterioration. Do not use if device appears damaged.
- Exercise care during operation to prevent any liquids or objects from spilling or falling into the product's housing.
- Always keep the alloy tank covered to minimize dross formation.
- After initial set-up keep alloy melter power on continuously.
- Do not allow alloy level to drop below 25mm (1") from bottom of alloy tank.
- Do not allow alloy level to exceed 50mm (2") from top of alloy tank.

1. Place alloy into melter with tongs one ingot at a time until appropriate level is achieved.

**⚠ WARNING**

- Do not add any other material, nor mix alloy types in the alloy melter tank.
- Hot alloy poses a burn hazard, do not touch.
- To remove alloy dross (oxidized alloy): Use ladle to extract and dispose in appropriate recycling bin.

2. Allow alloy tank temperature to stabilize within set range.
  - 79-82°C (175-180°F) for low temperature alloy.
  - 104-107°C (220-225°F) for medium temperature alloy.
3. Clear valve dispenser:
  - Place heat resistant vessel directly under nozzle.
  - Slowly open valve until alloy flows freely.
  - If alloy does not flow freely, close valve.
  - Increase valve temperature by one.
  - Wait 15 minutes and repeat.
4. Place Flat-Foam™ mold under nozzle and fill as appropriate.
5. Allow the alloy cast to cool to room temperature prior to handling.

**REPROCESSING****⚠ WARNING**

- Users of this product have an obligation and responsibility to provide the highest degree of safety hazard control to co-workers and themselves. To avoid safety hazards, wear all recommended protective equipment and follow safety control policies established by your facility.

1. Clean loose particles with a HEPA vacuum.
2. Clean with lead and heavy metal equipment cleanser or towelettes.
3. Wipe dry with lint-free cloth or allow to air dry prior to use.

**MAINTENANCE**

NOTE: • Inspect device prior to use for signs of damage and general wear.  
• If device does not function per instructions, contact CIVCO.

**DISPOSAL****⚠ WARNING**

- Ensure device is free of alloy before disposal.
- Contact CIVCO to recycle per EU Directives.

**STORAGE**

NOTE: For long term storage: Empty alloy from melter.

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COPYRIGHT © 2020 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
- AB Direktifleri uyarınca geri dönüřtürmek için CIVCO ile iletiřime gein.