Immobilizers

Dual Leg Immobilizer | GU Immobilizer

INSTRUCTIONS FOR USE (IFU)

CQ MEDICALTM

CE

INTENDED USE

The device is intended to aid in patient positioning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.


⚠ WARNING

- Do not use if device appears damaged.
- Do not reposition device with patient on it.
- Do not use Lok-Bar™ to secure positioning devices.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Ensure patient remains still for duration of setup and treatment.
- Product shape may have slight variations; the same immobilizer must be used for a patient throughout setup and treatment cycles.
- Ensure device is secure prior to use.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

DUAL LEG IMMOBILIZER

1. Place immobilizer onto Lok-Bar™ pins (Refer to Lok-Bar™ instruction) at appropriate position (A, B or C).
2. Position patient as dictated by setup protocol.
3. Secure patient feet with strap (provided).


- I. Immobilizer
- II. Lok-Bar™
- III. Lok-Bar™ Pins

GU IMMOBILIZER

1. Place immobilizer base onto Lok-Bar™ pins (Refer to Lok-Bar™ instruction) at appropriate location.
2. Release immobilizer pins to adjust cushions to appropriate position.

3. Position patient as dictated by setup protocol.

REPROCESSING

WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Clean device with a mild, non-abrasive cleaner as needed.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. LOK-BAR IS A TRADEMARK OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
**MEDTEC LLC**