# CT Simulation Laser QA Device

---

REFERENCE GUIDE

**CIVCO**  
Radiotherapy™

**INTENDED USE**

The device is intended to provide secondary radiologic quality assurance of CT simulators and room lasers following completion of primary system calibration.

**CAUTION**

Federal (United States) law restricts this device to sale by or on the order of a physician.


**⚠ WARNING**

- Do not use if device appears damaged or base is warped.
- Do not use Lok-Bar™ to secure devices.

- NOTE:
- Distance between side pegs from each center peg hole is 250mm +/- 1mm.
  - Verify 125mm +/- 0.5mm distance from center peg hole of middle peg to center peg hole of side pegs.

**CT SIMULATION LASER QA DEVICE SETUP**

1. Place device onto Lok-Bar™ pins.

**USING CT SIMULATION LASER QA DEVICE****IDENTIFY SCAN PLANE WITHIN GANTRY OPENING**

1. Adjust table height to align horizontal side gantry and vertical side-wall lasers with horizontal peg holes.
2. Adjust tabletop to align overhead axial gantry laser with vertical peg holes, side gantry lasers with horizontal peg holes, and vertical side lasers with horizontal peg hole.
3. Acquire single axial scan with 1 to 2 mm slice to ensure properly aligned device.

NOTE: Image should clearly reveal cross in each peg.

**ENSURE PARALLEL AND ORTHOGONAL GANTRY LASERS INTERSECT IN CENTER OF SCAN PLANE**

1. Raise and lower tabletop to ensure vertical side gantry lasers and overhead axial gantry laser always intersect side peg holes.
2. Raise and lower tabletop to ensure overhead sagittal gantry laser intersects center peg hole of middle peg.
3. Move tabletop in and out of gantry to ensure horizontal side gantry lasers always intersect side peg holes.
4. Move tabletop in and out of gantry to ensure overhead sagittal gantry laser intersects center peg hole of middle peg.

**ENSURE ACCURATELY SPACED VERTICAL SIDE-WALL LASERS**

1. Retract tabletop away from gantry to predefined separation between gantry and vertical side-wall lasers using digital longitudinal indicator.
2. Ensure vertical side-wall lasers bisect side peg holes.

**ENSURE PARALLEL AND ORTHOGONAL VERTICAL SIDE-WALL LASERS**

1. Raise and lower tabletop to ensure vertical side-wall lasers always intersect side peg holes.
2. Move tabletop in and out of gantry to ensure vertical side-wall lasers intersect side peg holes.

**ENSURE ORTHOGONAL OVERHEAD SAGITTAL LASERS**

1. Move tabletop in and out of gantry to ensure full overhead sagittal laser beam intersects center peg hole of middle peg.

**ENSURE ACCURATE, LINEAR AND REPRODUCIBLE OVERHEAD SAGITTAL LASER MOVEMENT**

1. Place ruler (provided) against 3 pegs and align overhead sagittal laser with mid mark of ruler.
2. Ensure overhead sagittal laser position reads zero.
3. Move laser side to side to verify laser motion accuracy.

**ENSURE LEVEL AND ORTHOGONAL TABLETOP**

1. Index device to most superior notches closest to gantry.
2. Align with gantry lasers and acquire a single slice through device.
3. Index device to most inferior notches within scan range.
4. Align with gantry lasers and acquire a single slice through device.
5. Measure location of center hole on both images.

NOTE: Measurements should be within 2mm agreement. Deviations on x and y axis greater than 2mm will indicate:

X axis: Tabletop is twisted or rotated.

y axis: Longitudinal roll in tabletop.

**ENSURE ACCURATE GANTRY MOVEMENT**

1. Tilt gantry and then return to vertical position.

NOTE: Alignment of vertical side gantry lasers should remain within 1mm from side peg holes.

**SCAN LOCALIZATION FROM SCOUT IMAGE**

1. Acquire horizontal scout image of device.
2. Use scout image to program scanner to acquire image through center of vertical holes.

NOTE: Image should clearly reveal holes in pegs.

REPROCESSING **WARNING**

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*

1. Remove visible contaminants using a soft damp cloth.
2. Wipe dry with lint-free cloth or allow to air dry.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2023 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. CIVCO RADIOTHERAPY AND CIVCO RT ARE TRADEMARKS OF MEDTEC LLC. LOK-BAR IS A TRADEMARK OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.