IPPS™ Overlay

An indexed patient positioning system for use with CTI PET/CT System.

REFERENCE GUIDE

CIVCO
Radiotherapy™

CE

INTENDED USE

The device is intended to aid in patient positioning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

- Do not use if device appears damaged.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Ensure device is secure prior to use.


NOTE:

- If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.
- Overlay weight: 30 lbs (14 kg). For total table load rating, including accessories, refer to your system's manual.

ATTACHING DOCKING STATION TO CRADLE

NOTE: MT-IL-3007 Tools Required: 3/8" hex wrench and level (*not provided*).


1. Remove all accessories from cradle.

- I. Cradle
- II. Docking Station
- III. Overlay
- IV. Side Support Blocks
- V. Assembly Latch
- VI. Shim

2. Fasten docking station and 1/8" shim loosely on the foot end of the table using screws and washers (*provided*).

ATTACHING OVERLAY

1. Place overlay onto cradle with side support blocks over anterior edge.

NOTE: May need to loosen side support blocks to fit around cradle edges.

2. Insert overlay into docking station.
3. Clip assembly latch into overlay and lock into place.
4. Ensure overlay is level to within 1° laterally and longitudinally.
5. Tighten all screws once positioning is verified.

NOTE:

- If vertical height adjustment is necessary; remove latch and adjust by adding or removing shims.
- If latch adjustment is necessary, loosen docking station screws, move longitudinally, retighten screws.


USING Lok-Bar™

WARNING

- Do not use Lok-Bar™ to secure positioning devices.


1. Snap Lok-Bar™ at appropriate longitudinal position. Ensure Lok-Bar™ is secure.

2. Place positioning accessory onto Lok-Bar™ pins.

- I. Overlay
- II. Lok-Bar™ Pins
- III. Lok-Bar™ Pins
- IV. Interloc® notches
- V. Scan Line

NOTE: Spacing between divergent wires (2) can be used as a reference for slice location when viewing multiple images containing axial slices.


- I. 24GA divergent wire
- II. 5° angle

NOTE: For optimal imaging, begin scanning after scan line.

⚠ WARNING

- Avoid pinch points when overlay is in motion by ensuring patient does not grasp overlay sides.
- Potential pinch point occurs between overlay and gantry when scanning with the table in lowest position.
- Potential pinch point occurs between overlay and table when retracting overlay.

- I. Gantry
- II. Overlay
- III. Table

REPROCESSING

⚠ WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2022 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. CIVCO RADIOTHERAPY AND CIVCO RT ARE TRADEMARKS OF MEDTEC LLC. IPPS AND LOK-BAR ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.