Shoulder Retractors

INSTRUCTIONS FOR USE (IFU)

INTENDED USE

The device is intended to aid in patient positioning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.


⚠ WARNING

- Do not use if device appears damaged.
- Do not reposition device with patient on it.
- Do not use excessive force on device. Injury or device damage may occur.
- Ensure device/patient does not move for duration of set-up and treatment.
- Device is for use with CQ Medical specified accessories only.
- Ensure device is secure prior to use.
- Verify patient position with completed setup sheet prior to treatment.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at www.CQmedical.com.
- Verification image should be taken to confirm patient position prior to treatment.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

USING SHOULDER RETRACTORS

1. Insert rope/handle assemblies through sockets on foot base and lock indexing balls into slots.

2. Position patient's feet against base with patient lying down.
3. Assist patient with handle bars or wrist loops.
4. Pull ropes to retract patient's shoulders. Lock indexing balls into place and record locations.

⚠ WARNING

- Switching strap styles will result in different lengths.

REPROCESSING**⚠ WARNING**

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Clean non-fabric components with mild, non-abrasive cleanser. Replace fabric wrist loops if they become contaminated. Optional blue nylon wrist loop may be cleaned by wiping all surfaces with common germicidal or antiseptic, such as alcohol.

NOTE: If device becomes contaminated or can no longer be cleaned, discontinue use and replace.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
**MEDTEC LLC**