Shoulder Retractors

INSTRUCTIONS FOR USE (IFU)

INTENDED USE

The device is intended to aid in patient positioning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.


WARNING

- Do not use if device appears damaged.
- Do not reposition device with patient on it.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at www.CQmedical.com.
- Verify patient position with completed setup sheet prior to treatment.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Ensure device/patient does not move for duration of set-up and treatment.
- Verification image should be taken to confirm patient position prior to treatment.
- Ensure device is secure prior to use.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

USING SHOULDER RETRACTORS

1. Position patient's feet into Feetfix™ Couch Support and cushion as needed for desired comfort.

2. Insert fixing knob into appropriate indexing hole on Feetfix™ Couch Support.
3. Assist patient with wrist loops.

NOTE: Ensure fixing knob is pushed far enough into the holes to prevent it from pulling out.

REPROCESSING

WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol. Connecting strap may be laundered as needed. If device is contaminated and can no longer be cleaned, discontinue use and replace.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. FEETFIX IS A TRADEMARK OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MEDTEC LLC