


Posifix®/Type-S™ Overlay

Head-neck-shoulders patient immobilization system for use with I-Beam

REFERENCE GUIDE

CIVCO

CE

- I. Couchtop
- II. Posicast® Overlay
- III. Locking mechanisms (2)
- IV. Type-S™ Overlay
- V. Overlay pins

INTENDED USE

Device is intended for patient positioning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Do not use if device appears damaged.
- Do not allow patient to reposition themselves.
- Do not reposition device with patient on it.
- Do not exceed 55 lbs (25 kg) extended past end of table.
- This product is manufactured with carbon fiber, and rough handling may cause damage resulting in fibers coming off. If this happens, stop using the product to prevent fibers from getting in contact with the patient's or user's skin.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

USING OVERLAY

1. Engage overlay pins into indexing holes on couchtop.
2. Ensure overlay is secured by locking mechanisms.
3. Place desired headrest (*refer to headrest instruction*).
4. Position patient as dictated by setup protocol. It is recommended assistance be given in slowly positioning patient downward (*hips to head*) onto overlay to ensure maximum working load of 55 lbs (25 kg) is not exceeded.

USING THERMOPLASTICS WITH TYPE-S™ OVERLAY

1. Place shim (*provided*) onto baseplate pins, if desired.
2. Form and attach thermoplastic mask (*refer to Thermoplastic instruction*).

USING POSICAST® THERMOPLASTICS WITH POSICAST OVERLAY

1. Form and attach thermoplastic mask (*refer to Thermoplastic instruction*).

REPROCESSING**⚠ WARNING**




- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Disinfect with a mild, non- abrasive cleanser.
2. Wipe dry with lint- free cloth or allow to air dry prior to reassembly and use.

MAINTENANCE**⚠ WARNING**

- Inspect locking mechanisms prior to each use for signs of damage.

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2020 ALL RIGHTS RESERVED. CIVCO AND POSICAST ARE REGISTERED TRADEMARKS OF CIVCO MEDICAL SOLUTIONS. TYPE-S IS A TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.