Posifix® Spacers & Wedges

INSTRUCTIONS FOR USE (IFU)

INTENDED USE

107101, 107102, 107103, 107104: Fixation and repositioning of the head and neck area of the patient during radiotherapy and diagnostics.

107101MR, 107102MR, 107103MR, 107104MR: The device is indicated to aid in supporting, positioning, and/or immobilization of adult and pediatric patients undergoing radiation therapy of the head, brain, neck, and spine including radiosurgery and electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

- Do not use if device appears damaged.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- To avoid patient setup error, combined devices must be stacked in same order throughout patient treatment cycle. Record how individual devices are stacked.
- Ensure device/patient does not move for duration of set-up and treatment.
- Verification image should be taken to confirm patient position prior to treatment.
- Ensure device is secure prior to use.
- Do not reposition device with patient on it.
- Devices must be combined in the same order for each treatment.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

MRI SAFETY INFORMATION

MR Safe

Devices are MR Safe:

- 107101MR
- 107102MR
- 107103MR
- 107104MR
- 107106MR

USING SPACER & WEDGE

1. Place desired spacers and/or wedges with optional Posifix® head and neck treatment device.
2. Place Posifix® headrest into rectangular indentation on desired spacer or wedge, compressing foam headrest inwards on the sides as needed to ensure headrest is fully seated.

NOTE: Use the most superior thermoplastic docking location.

REPROCESSING

WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. POSIFIX IS A REGISTERED TRADEMARK OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MEDTEC LLC