Body Pro-Lok™ Platform

INSTRUCTIONS FOR USE (IFU)

CE


Towards Gantry

Towards Gantry
朝向扫描架

Prema kućištu
Směrem k nosné konstrukci
Mod ophængning
Naar rijbrug
Pöytää kohti
Couvercle de portique
Zum Gestell hin
Προς πλαίσιο
A gantry irányába
Verso il gantry
構台方向

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*Images are not to scale and are intended for display purposes only.

INTENDED USE

The device is intended to aid in patient positioning and immobilizing during radiotherapy procedures.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING


- Do not use if device appears damaged.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Ensure all Body Pro-Lok™ accessories are properly attached using T-pin Lok-Bars or bridge latching mechanisms.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at www.CQmedical.com.
- Verify patient position with completed setup sheet prior to treatment.
- Ensure patient remains still for duration of setup and treatment.
- Monitor patient constantly when using accessories to ensure patient comfort.
- Ensure clearance is maintained between Body Pro-Lok™ system and all linac or imaging equipment before using system.
- Do not attempt to move or carry platform by accessory bridges.
- Ensure hands, feet and other extremities are kept away from moving parts.
- Do not use Lok-Bar™ to secure devices.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

BODY PRO-LOK™ PLATFORM**⚠ WARNING**

- Do not exceed patient weight of 500 lbs (226.8 kg) while stationary.
- Ensure Body Pro-Lok™ Platform is securely resting on tabletop using latching transfer system or Lok-Bar™.
- Ensure Body Pro-Lok™ Platform is level to prevent patient sliding or falling.
- For more than one platform, ensure all platform lateral scale indexing markers align. If platforms have opposite orientation, contact your CQ Medical representative.

1. Place Body Pro-Lok™ Platform onto Lok-Bar™ pins.

- I. Body Pro-Lok™ Platform
- II. Lok-Bar™ (2)
- III. Lok-Bar™ Pins


REPROCESSING**⚠ WARNING**

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. BODY PRO-LOK, CLAM-LOK AND LOK-BAR ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MEDTEC LLC