# MRI Type-S™

Head-neck-shoulders patient immobilization system

---

INSTRUCTIONS FOR USE (IFU)

**CQ MEDICAL™**

CE

**INTENDED USE**

The device is indicated to aid in supporting, positioning, and/or immobilization of adult and pediatric patients undergoing radiation therapy of the head, brain, neck, and spine including radiosurgery and electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.

**CAUTION**

Federal (United States) law restricts this device to sale by or on the order of a physician.

**⚠ WARNING**

- Verify all angles of treatment, attenuation characteristics and WET values prior to treating patients.
- Do not use if device appears damaged.
- Do not use Lok-Bar™ to secure positioning devices.
- Do not allow patient to reposition themselves, such as arching back.
- Do not reposition device with patient on it.
- Ensure device is secure prior to use.
- MTMRHNTYPES: Do not apply a load exceeding 55 lbs (24.9 kg).

**NOTE:** If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

**MRI SAFETY INFORMATION**

MR Safe

- Device is MR safe.


**USING TYPE-S™ OVERLAY**

1. Place Type-S™ overlay onto Lok-Bar™ pins (refer to Lok-Bar™ instruction).

**⚠ WARNING**

- Do not extend inscribed line past end of table.

- NOTE:**
- Velcro strap is not required when device is not extended.
  - Attach overlay to table using Velcro strap before positioning patient.
2. Wrap velcro strap around Type-S™ overlay and table.
  3. Place desired headrest onto Type-S™ (refer to headrest instruction).
  4. Position patient as dictated by setup protocol.

**USING THERMOPLASTICS**

- NOTE:**
- Devices are compatible with CQ Medical Type-S™ thermoplastics (refer to Thermoplastic instruction).
  - Only use shim when first forming thermoplastic mask.

1. Place shim (provided) onto baseplate pins, if desired.
2. Form and attach thermoplastic mask (refer to Thermoplastic instruction).

**REPROCESSING****WARNING**

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*

1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol.

**MAINTENANCE**

**NOTE:** Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. LOK-BAR, TYPE-S AND VAC-LOK ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
**MEDTEC LLC**