Universal Couchtop

for use with Elekta Precise, Elekta Synergy, Siemens TXT, Siemens ZXT, Varian Exact and CQ Medical Protura systems

INSTRUCTIONS FOR USE (IFU)

INTENDED USE

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

ATTACHING UNIVERSAL COUCHTOP™ TO ELEKTA PRECISE AND SYNERGY

NOTE: Tools Required: 6mm Hex Wrench, Level (*not provided*).

1. Remove center pivot and latch.

2. Place couchtop on table and align attachment holes (6).
3. Attach couchtop to table using hardware (provided). Torque bolts to 21 N·m (15.5 lb·ft).

4. Ensure couchtop is level to within 0.25° laterally and longitudinally.
- NOTE: If couchtop is not level, remove couchtop and shim accordingly until level is achieved.
5. Snap Lok-Bar™ at appropriate longitudinal position. Ensure Lok-Bar™ is secure.
6. Attach positioning device onto Lok-Bar™ pins.

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⚠ WARNING


- Do not use if device appears damaged.
- Verify all angles of treatment, attenuation characteristics and WET values prior to treating patients. Refer to Technical Data Sheet at www.CQmedical.com for WET values.
- Do not exceed a patient weight of 440 lbs (200 kg).

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
ATTACHING UNIVERSAL COUCHTOP™ TO SIEMENS TXT AND CQ MEDICAL PROTURA™

- NOTE:
- Tools Required: 8mm Hex Wrench, screwdriver and leveling tool (*not provided*).
 - Remove existing top.
 - Use hardware provided with Protura™ for installation on Protura™ unit.

1. Place spacers on pedestal at attachment locations.

2. Attach couchtop to table using hardware (*provided*).

3. Ensure couchtop is level to within 0.25° laterally and longitudinally.


NOTE: If couchtop is not level, remove couchtop and shim accordingly until level is achieved.

ATTACHING UNIVERSAL COUCHTOP™ TO SIEMENS ZXT

- NOTE:
- Tools Required: 5mm and 8mm Hex Wrenches, screwdriver and leveling tool (*not provided*).
 - Remove head extension, half panels, grid, center panel and all other accessories.

1. Rotate table 45° to expose bolts and lift steel frame off of table.
2. Remove center pivot bearing from frame and refasten on interface plate.
3. Attach interface plate to table using hardware (*provided*).
4. Attach interface plate to couchtop using hardware (*provided*).
5. Rotate and lock couchtop onto table using latch handle.
6. Ensure couchtop is level to within 0.25° laterally and longitudinally.

NOTE: If couchtop is not level, remove couchtop and shim accordingly until level is achieved.

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- Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol.

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
ATTACHING UNIVERSAL COUCHTOP™ TO VARIAN EXACT


NOTE: Tools Required: 5/16", 4mm and 8mm Hex Wrenches and Level (not provided).

1. Remove existing table and components, including frame.
2. Attach rectangular spacer to interface plate using hardware (provided) and place round spacers onto pedestal.


- I. Pedestal
- II. Interface Plate
- III. Round Spacers
- IV. Rectangular Spacer
- V. M5 Screws

3. Attach interface plate to pedestal using hardware (provided).

- I. 3/8" Screws

4. Attach couchtop to interface using hardware (*provided*).

5. Ensure couchtop is level to within 0.25° laterally and longitudinally.

NOTE: If couchtop is not level, remove couchtop and interface plate; shim accordingly until level is achieved.

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INTENDED USE

- MTIL6720: The device is indicated to aid in supporting and positioning adult and pediatric patients undergoing radiation therapy including electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.
- MTIL6020, MTIL6021, MTIL6121, MTIL6220, MTIL6601, MTIL6604, MTIL6605, MTIL6606, MTIL6607, MTIL6608, MTIL6610: The device is intended to aid in patient positioning.

CAUTION


Federal (United States) law restricts this device to sale by or on the order of a physician.

USING EXTENSIONS (IF APPLICABLE)**⚠ WARNING**


- Do not use if device appears damaged.
- Verify all angles of treatment, attenuation characteristics and WET values prior to treating patients. Refer to Technical Data Sheet at www.CQmedical.com for WET values.
- It is not recommended to treat using posterior oblique beams in the lateral area of the couch/extension junction, where the lock mechanism is contained.
- Do not allow patient to reposition themselves, such as arching back.
- Type-S™ Extension: Do not apply a load exceeding 150 lbs (68 kg).
- Rectangular Extension: Do not apply a load exceeding 150 lbs (68 kg).
- Posifix® Extension: Do not apply a load exceeding 150 lbs (68 kg).
- Palliative Extensions: Do not apply a load exceeding 150 lbs (68 kg).
- Ensure extension is in locked position by lifting gently on superior end.
- Avoid pinch points when overlay is in motion by ensuring patient does not grasp overlay sides.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

1. Place extension onto superior end of table, ensuring extension is secure.
2. Place desired positioning device onto extension (refer to appropriate positioning device manual).
3. Position patient as dictated by setup protocol.

**REMOVING EXTENSIONS (IF APPLICABLE)**


1. Rotate lever and lift extension upward to remove extension.

USING VARIAN INTERFACE EXTENSION (IF APPLICABLE) **WARNING**

- Do not allow patient to reposition themselves, such as arching back.
- Inspect device to ensure support braces are not loose, cracked or have become detached.
- Do not reposition device with patient on it.
- Do not apply a load exceeding 55 lbs (24.9 kg).

1. Place Varian interface extension onto superior end of table, ensuring extension is secure.
2. Latch extension onto Varian interface, ensuring extension is secure.
3. Place desired positioning device onto extension (refer to appropriate positioning device manual).
4. Position patient as dictated by setup protocol.

REPROCESSING **WARNING**

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. IPPS, LOK-BAR, PRODIGY, AND TYPE-S ARE TRADEMARKS OF MEDTEC LLC. POSIFIX IS A REGISTERED TRADEMARK OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MEDTEC LLC