Dual Function Vac-Lok™ Vacuum Compressor

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| Test / Standard | Emissions Class and Group / Immunity Test Level |
|---|--|
| CISPR 11, Conducted Emissions | Class B; Group 1 |
| CISPR 11, Radiated Emissions | Class B; Group 1 |
| IEC 61000-3-2, Harmonic Current Emissions | Class A |
| IEC 61000-3-3, Voltage Changes, Voltage Fluctuations, and Flicker | All Parameters But Long-Term Flicker |
| IEC 61000-4-2, Electrostatic Discharge Immunity | ±8kV contact ±2kV, ±4kV, ±8kV, ±15kV air |
| IEC 61000-4-3, Radiated RF Immunity | See Note Below |
| IEC 61000-4-3, Proximity fields from RF wireless equipment | |
| IEC 61000-4-4, Electrical Fast Transient / Burst Immunity | |
| IEC 61000-4-5, Surge Immunity | ±0.5kV, ±1kV Line-to-Line ±0.5kV, ±1kV, ±2kV Line-to-Ground |
| IEC 61000-4-6, Immunity to Conducted Disturbances, Induced by Radio Frequency Fields | See Note Below |
| IEC 61000-4-8, Magnetic Field Immunity | |
| IEC 61000-4-11, Voltage Dips | |
| IEC 61000-4-11, Voltage Interruptions | |
| Note: The EUT contains no active electrical or electronic components to malfunction from the immunity tests listed above. | |

| Replacement Parts - Contact CQ Medical for service. | |
|---|-----------------------------|
| MTVP1 | Power Cord EJCFPCSETK3 |
| MTVPExA | Power Cord EJCFPCSETA |
| MTVPExC | Power Cord EJCFPCSETC |
| MTVPExD | Power Cord EJCFPCSETD |
| MTVPExG | Power Cord EJCFPCSETG |
| MTVPExH | Power Cord EJCFPCSETH |
| MTVPFT | Vacuum Pump Foot Pedal |
| MTVPHD | Vacuum Pump Hand Control |
| MTVPHS | Vacuum Hose |

PRIOR TO USE:

- The user is responsible for ensuring the application and use of the device does not compromise the patient contact rating of any equipment used in the vicinity of, or in conjunction with, the system.
- The use of accessory equipment and/or hardware not complying with the equivalent product safety and EMC requirements of this product may lead to a reduced level of safety and/or EMC performance of the resulting system. Consideration relating to the choice of accessory equipment used with this product shall include:
 - The use of the accessory in the patient vicinity.
 - Evidence the safety certification of the accessory has been performed in accordance with the appropriate IEC 60601-1 Harmonized National Standards.
 - Evidence the EMC certification of the accessory has been performed in accordance to the IEC 60601-1-2 Harmonized National Standards.
- Observe all safety precautions recommended by the accessory equipment manufacturer in the user documentation provided with the equipment.
- Keep this manual with device for reference as required.

INTENDED USE

Device is intended to inflate and/or deflate Vac-Lok™ cushions which is intended for patient positioning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.


⚠ WARNING

- Routinely check vacuum compressors and power cords for signs of deterioration.
- Do not use if device appears damaged.
- Do not cover ventilation openings.
- Do not attempt to service device.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Vac-Lok™ Vacuum Compressor is for use with CQ Medical Vac-Lok™ Cushions only.
- Do not continuously inflate or deflate Vac-Lok™ cushion for more than 2 minutes per hour.
- Equipment is not to be modified without CQ Medical authorization.
- To avoid risk of electric shock, only connect equipment to a supply mains with protective earth.
- Users of device should ensure they have read and understand instructions for use prior to operating device.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Dual Function Vac-Lok™ Vacuum Compressor, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

- NOTE:
- If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.
 - Place female valve connector onto cushion valve to let air into Vac-Lok™ cushion without using pump.
 - Dual Action Vac-Lok™ Compressor is suitable for use in Radiation Oncology planning and treatment environments.

USING DUAL FUNCTION VAC-LOK™ VACUUM COMPRESSOR

1. Verify Vac-Lok™ Vacuum compressor's ON/OFF switch is OFF.
2. Plug power cord into Vac-Lok™ Vacuum Compressor.
3. Plug hand or foot pendant into Vac-Lok™ Vacuum Compressor.
4. Plug power cord into appropriate GFCI (ground fault circuit interrupter) receptacle.
5. Plug vacuum hose into Vac-Lok™ vacuum compressor.
6. Turn device ON.
7. Attach vacuum hose to cushion valve (refer to Vac-Lok™ Cushion manual).

8. Briefly press INFLATE ↑ and DEFLATE ↓ buttons on hand or foot pendant to verify unit works in both modes (refer to Vac-Lok™ Cushion manual).

9. Deflate until audible alarm sounds.

NOTE: Dual Action Vac-Lok™ Compressor shall pull vacuum until the audible chime, resulting in a rigid Vac-Lok™ Bag. If performance is lost, due to electromagnetic disturbances, the Vac-Lok™ Bag will not be rigid.

10. Verify device's ON/OFF switch is OFF.
11. Disconnect vacuum hose from Vac-Lok™ cushion.

REPROCESSING

 **WARNING**

- *Allow device to dry before plugging into electrical socket.*

1. Clean device with a mild, non-abrasive cleaner as needed.

MAINTENANCE

 **WARNING**

- *Failure to use specified fuse could result in device damage or injury.*

NOTE: Inspect device prior to use for signs of damage and general wear.

DISPOSAL

 **WARNING**

- *Device should be disposed of in manner appropriate for electronic equipment.*

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| Test / Standard | Emissions Class and Group / Immunity Test Level |
|---|---|
| CISPR 11, Conducted Emissions | Class B; Group 1 |
| CISPR 11, Radiated Emissions | Class B; Group 1 |
| IEC 61000-3-2, Harmonic Current Emissions | Class A |
| IEC 61000-3-3, Voltage Changes, Voltage Fluctuations, and Flicker | All Parameters But Long-Term Flicker |
| IEC 61000-4-2, Electrostatic Discharge Immunity | ±8kV contact ±2kV, ±4kV, ±8kV, ±15kV air |
| IEC 61000-4-3, Radiated RF Immunity | See Note Below |
| IEC 61000-4-3, Proximity fields from RF wireless equipment | |
| IEC 61000-4-4, Electrical Fast Transient / Burst Immunity | |
| IEC 61000-4-5, Surge Immunity | |
| IEC 61000-4-6, Immunity to Conducted Disturbances, Induced by Radio Frequency Fields | See Note Below |
| IEC 61000-4-8, Magnetic Field Immunity | |
| IEC 61000-4-11, Voltage Dips | |
| IEC 61000-4-11, Voltage Interruptions | |
| Note: The EUT contains no active electrical or electronic components to malfunction from the immunity tests listed above. | |

| Reservdelar - Kontakta CQ Medical för service. | |
|--|--------------------------------|
| MTVP1 | Nätsladd EJCFPCSETK3 |
| MTVPExA | Nätsladd EJCFPCSETA |
| MTVPExC | Nätsladd EJCFPCSETC |
| MTVPExD | Nätsladd EJCFPCSETD |
| MTVPExG | Nätsladd EJCFPCSETG |
| MTVPExH | Nätsladd EJCFPCSETH |
| MTVPFT | Fotpedal till vakuumpump |
| MTVPHD | Handreglage till vakuumpump |
| MTVPHS | Vakuumslang |

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. VAC-LOK IS A TRADEMARK OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
**MEDTEC LLC**