Posifix® Baseplate

INSTRUCTIONS FOR USE (IFU)

INTENDED USE

The device is intended to aid in patient positioning during radiotherapy and diagnostics treatments.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

- Do not use if device appears damaged.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Do not reposition device with patient on it.
- Do not use Lok-Bar™ to secure positioning devices.
- Only use compatible CQ Medical accessories in combination with this device.
- Ensure device is secure prior to use.

- NOTE:
- If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.
 - Baseplates require the use of a headrest (*refer to Headrest instruction*).

USING ACRYLIC POSIFIX® BASEPLATE

MRI SAFETY INFORMATION

MR Safe

- Acrylic Posifix® Baseplate is safe for use within MR environment.

1. Place Posifix® baseplate onto Lok-Bar™ pins (*refer to Lok-Bar™ instruction*).
2. Place desired headrest onto baseplate.
3. Position patient onto baseplate for thermoplastic setup (*refer to thermoplastic instruction*).
4. Attach quickfasteners to baseplate to secure thermoplastic.

USING CARBON FIBER POSIFIX® BASEPLATE

MRI SAFETY INFORMATION

MR UNSAFE

- Carbon Fiber Posifix® Baseplate is not safe for use within MR environment.

1. Place Posifix® baseplate onto Lok-Bar™ pins (*refer to Lok-Bar™ instruction*).

NOTE: Baseplate fixed with aluminium rail will differ in position of patient in longitudinal direction by 44mm from a baseplate fixed with a Lok-Bar™.

2. Place desired headrest onto baseplate.
3. Position patient onto baseplate for thermoplastic setup (*refer to thermoplastic instruction*).
4. Attach quickfasteners to baseplate to secure thermoplastic.

REPROCESSING

WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Clean device with a mild, non-abrasive cleaner as needed.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. POSIFIX IS A REGISTERED TRADEMARK OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
**MEDTEC LLC**