# New Horizon™ Prone Breastboard


---

REFERENCE GUIDE

**CIVCO**  
Radiotherapy™

CE

- |   |   |
|---|---|
| <ul style="list-style-type: none"> <li>I. Overlay/Couchtop</li> <li>II. Lok-Bar™</li> <li>III. Spacer Assembly</li> <li>IV. Inferior Section</li> <li>V. Superior Section</li> <li>VI. Center Support Section</li> <li>VII. Docking Button</li> </ul> | <ul style="list-style-type: none"> <li>VIII. Breast Wedge</li> <li>IX. Wedge Slide Tray</li> <li>X. Ipsilateral Insert</li> <li>XI. Head Cushion</li> <li>XII. Handle</li> <li>XIII. Bridge Rail</li> <li>XIV. Cushion</li> </ul> |
|---|---|

**INTENDED USE**

The device is intended to aid in patient positioning.

**CAUTION**

Federal (United States) law restricts this device to sale by or on the order of a physician.

**⚠ WARNING**

- Do not use if device appears damaged.
- Do not allow patient to reposition themselves.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Do not use Lok-Bar™ to secure positioning devices.
- Ensure device is secure prior to use.
- Only use CIVCO Radiotherapy compatible accessories.
- Ensure all accessories are attached securely prior to treatment.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at [WWW.CIVCORT.COM](http://WWW.CIVCORT.COM).
- Verify patient position with completed setup sheet prior to treatment.
- Do not exceed patient weight of 500 lbs (226 kg).
- Ensure patient remains still for duration of setup and treatment.
- Do not reposition device with patient on it.

- NOTE:
- If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.
  - Total Device weight: 60 lbs (27 kg)
    - Superior section: 27 lbs (12 kg)
    - Inferior section: 28 lbs (13 kg)
    - Center section: 5 lbs (2 kg)

**USING NEW HORIZON™ PRONE BREASTBOARD**

1. Attach spacer assemblies equally to bottom of inferior and superior sections as needed.
2. Place inferior and superior sections onto Lok-Bar™ pins.
3. Slide center support section onto docking buttons between inferior and superior sections with opening on desired side for treatment.
4. Position desired breast wedge cushion and wedge slide tray (*if desired*) on contralateral breast side.
 

NOTE: Adjust wedge slide tray during patient positioning to move excess tissue laterally, away from the treatment field.
5. Place ipsilateral insert into treatment side of center support (*if desired*) and cover with corresponding size cushion.
6. Position head cushion on superior section, record indexing on setup sheet.
 

NOTE: Head cushion can be moved superior or inferior as well as anterior or posterior as desired for patient comfort.
7. Position handles in a way that is comfortable for the patient, record indexing on setup sheet.
8. Attach SBRT Body Pro-Lok™ Bridge (*optional*) to bridge rail if desired (*refer to SBRT Body Pro-Lok™ instruction*).
9. Position the patient as dictated by setup protocol.





**REPROCESSING****⚠ WARNING**

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol.

**MAINTENANCE**

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2022 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. CIVCO RADIOTHERAPY AND CIVCO RT ARE TRADEMARKS OF MEDTEC LLC. BODY PRO-LOK, HORIZON AND LOK-BAR ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.