Universal Couchtop Extensions

for use with UCT LE

INSTRUCTIONS FOR USE (IFU)

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INTENDED USE

The device is intended to support and aid in positioning a patient during radiologic and other medical procedures.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

- Do not use if device appears damaged.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Do not apply a load exceeding 484 lbs (220 kg).
- Avoid pinch points when overlay is in motion by ensuring patient does not grasp overlay sides.

NOTE:
 • If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.
 • It is not recommended to treat using posterior oblique beams in the lateral area of the couch/extension junction, where the lock mechanism is contained.


USING EXTENSIONS


WARNING

- Type-S™ Extension: Do not apply a load exceeding 315 lbs (143.2 kg).
- Grid Extension/Rectangular Extension: Do not apply a load exceeding 392 lbs (178 kg).
- Posifix® Extension: Do not apply a load exceeding 315 lbs (143.2 kg).


1. Insert fingers of extension into receiver on pedestal side of couch.
2. Lower into place and ensure extension is locked in position by lifting gently on superior end.


3. Place desired positioning device onto extension (refer to appropriate positioning device manual).
4. Position patient as dictated by setup protocol.


NOTE: Line indicates where clear scanning/treating area ends due to increased density in structure of extension.

NOTE: Spacing between divergent wires (2) can be used as a reference for slice location when viewing multiple images containing axial slices.


NOTE: Patient isocenter must be located superior to painted line.

REMOVING EXTENSIONS

1. Press inferior end of release lever until it locks down.
2. Remove extension using both hands.


USING VARIAN INTERFACE EXTENSION (IF APPLICABLE)

WARNING

- Inspect device to ensure support braces are not loose, cracked or have become detached.
- Do not apply a load exceeding 315 lbs (143.2 kg) to the entire extension, including accessory.

1. Insert fingers of extension into receiver on pedestal side of couch.
2. Lower into place and ensure extension is locked in position by lifting gently on superior end.
3. Latch extension onto Varian interface, ensuring extension is secure.
4. Place desired positioning device onto extension (refer to appropriate positioning device manual).
5. Position patient as dictated by setup protocol.

REPROCESSING

WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. IPPS, LOK-BAR, PRODIGY, AND TYPE-S ARE TRADEMARKS OF MEDTEC LLC. POSIFIX IS A REGISTERED TRADEMARK OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MEDTEC LLC