Body Pro-Lok™ Rails-Only System

REFERENCE GUIDE

CIVCO

CE

 Towards Gantry	Towards Gantry 朝向扫描架 Prema kućištu Směrem k nosné konstrukci Mod ophængning Naar rijbrug Pöytää kohti Couvercle de portique Zum Gestell hin Προς πλαίσιο A gantry irányába Verso il gantry 構台方向 Gentrija virzienā Mot gantry W kierunku suwnicy Em direção à estrutura Em direcção à torre Spre suportul de scanare Protí podstavku Hacia la pasarela Mot gantry Kızağa Doğru
--	--

INTENDED USE

The device is intended to aid in patient positioning and immobilizing during radiotherapy procedures.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

**WARNING**


- Do not use if device appears damaged.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Ensure all Body Pro-Lok™ accessories are properly attached using T-pin Lok-Bars or bridge latching mechanisms.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at WWW.CIVCORT.COM.
- Verify patient position with completed setup sheet prior to treatment.
- Ensure patient remains still for duration of setup and treatment.
- Do not use Lok-Bar™ to secure devices.
- Monitor patient constantly when using accessories to ensure patient comfort.
- Ensure clearance is maintained between Body Pro-Lok™ system and all linac or imaging equipment before using system.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

RAILS- ONLY SYSTEM

1. Attach Rails- Only System to couchtop at appropriate longitudinal position. Ensure Rails- Only System is secure.

NOTE: • Note all index marks for desired position.
• Note direction of Rails- Only System, with respect to gantry.

**REPROCESSING****WARNING**

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2020 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. BODY PRO-LOK AND LOK-BAR ARE TRADEMARKS OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.