MR Posi-S™ Adapter

REFERENCE GUIDE

CIVCO
Radiotherapy™

CE

INTENDED USE

The device is intended to aid in patient positioning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Do not use if device appears damaged.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Do not reposition device with patient on it.
- Verification image should be taken to confirm patient position prior to treatment.
- Ensure device is secure prior to use.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.


MRI SAFETY INFORMATION

MR Safe


- MR Posi-S™ Adapter is MR Safe.
- In non-clinical testing, the image artifact caused by the device extends approximately 12.8mm from the edge of the MR Posi-S™ Adapter pocket when imaged with a gradient echo pulse sequence and a 3.0T MRI System.

USING MR POSI-S™ ADAPTER

1. Attach Type-S™ Baseplate onto tabletop. (Refer to Type-S™ instruction).
2. Place MR Posi-S™ Adapter onto Type-S™ pins.

3. Place desired Posifix® headrest onto baseplate.

I. Headrest

4. Position patient as dictated by setup protocol.

REPROCESSING

⚠ WARNING

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*

1. Clean product with mild, non-abrasive cleanser.

MAINTENANCE

NOTE: Device must be inspected before and after use for damage and wear.

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COPYRIGHT © 2022 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. CIVCO RADIOTHERAPY AND CIVCO RT ARE TRADEMARKS OF MEDTEC LLC. POSI-S, TYPE-S, AND POSIFIX ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.