# Fiducial Markers

---

REFERENCE GUIDE

**CIVCO**  
Radiotherapy™

CE  
1639

**INTENDED USE**

Device is intended to be placed into soft tissue or bone to establish a permanent, accurate, internal reference system.

**CAUTION**

Federal (United States) law restricts this device to sale by or on the order of a physician.

**⚠ WARNING**

- Do not use if device appears damaged.
- Before use, you should be trained in minimally invasive biopsy techniques.
- Patients with known sensitivity to the materials listed in the device material section (see below) may suffer an allergic reaction to this device. If there is a known sensitivity to a device material, it is the responsibility of the physician to evaluate the risk/benefit prior to the use of the device.
- Disposable components are packaged sterile and are single-use only.
- Product is sterilized via ethylene oxide. Do not use if integrity of packaging is violated or if expiration date has passed.
- Dispose of single-use components as infectious waste.
- Ensure fiducial markers are securely implanted.
- Ensure fiducial markers are as well separated as possible.
- Create an approximately equilateral triangle when implanting three or more fiducial markers.
- Ensure fiducial markers do not overlie each other in AP or lateral projection.
- Do not use in lung.
- Fiducial migration is possible after implantation; verify final marker placement 3-5 days after implantation.
- Do not reuse, reprocess or resterilize single use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.

**NOTE:** If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

**DEVICE MATERIALS**

PART NUMBER	MATERIAL
MTNW887805, MTNW887808, MTNW887809, MTNW887812, MTNW887854, MTNW887855	99.99% Gold
MTNW887825, MTNW887830, MTNW887835, MTNW887836, MTNW887837, MTNW887840, MTNW887851, MTNW887853, MTNW887860, MYNW887870, MTNW887871, MTNW887872, MTNW887873, MTNW887881, MTNW887882, MTNW887883	99.99% Gold, Bone Wax
MTNW887351, MTNW887358	Bone Wax

**USING BONE MARKERS**

**NOTE:** Bone Markers (MTNW887805) are not CE marked and not available in the European market.

1. Remove fiducial marker from package, making sure to use proper sterile technique.
2. Locate intended site using appropriate imaging technique.
3. Set drill depth.
4. Drill and place fiducial marker.
5. Backfill hole with bonewax.
6. Confirm final marker position(s) with desired imaging.

**USING SOFT TISSUE MARKERS**

1. Remove soft tissue marker from package, making sure to use proper sterile technique.
2. Locate intended site using appropriate imaging technique.
3. Direct needle to intended site.
4. Advance stylet to position marker at tip of needle, then retract needle over stylet to extrude marker.
5. Remove needle and repeat as necessary.
6. Confirm final marker position(s) with desired imaging.

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COPYRIGHT © 2023. ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. CIVCO RADIOTHERAPY AND CIVCO RT ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.