Body Pro-Lok™ Accessories

Body Pro-Lok™ System

INSTRUCTIONS FOR USE (IFU)

	Towards Gantry باتجاه القطار 朝向扫描架 Prema kućištu Směrem k nosné konstrukci Mod ophængning Naar rijbrug Portaali suunas Pöytää kohti Couvercle de portique Zum Gestell hin Προς πλαίσιο A gantry irányába Verso il gantry 構台方向 갠트리 방향 Gentrija virzienā Link rémo Mot gantry W kierunku suwnicy Em direção à estrutura Em direção à torre Spre suportul de scanare К гентрикК гентри Prema postolju Smerom k portálu Proti podstavku Hacia la pasarela Mot gantry ไปทางโครงสร้างลิ้นชัก Kızağa Doğru
---	--

INTENDED USE

The device is intended to aid in support and positioning of patients during MR, radiological and other procedures.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Do not use if device appears damaged.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Ensure all Body Pro-Lok™ accessories are properly attached using Lok-Bars or bridge latching mechanisms.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at www.CQmedical.com.
- Verify patient position with completed setup sheet prior to treatment.
- Ensure patient remains still for duration of setup and treatment.
- Do not use Lok-Bar™ to secure devices.
- Monitor patient constantly when using accessories to ensure patient comfort.
- Ensure clearance is maintained between Body Pro-Lok™ system and all linac or imaging equipment before using system.
- Do not attempt to move or carry platform by accessory bridges.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

MRI SAFETY INFORMATION

MR Safe

- Forehead Brace is MR Safe.
- Shoulder Bridge 2 is MR Safe.
- Rails-Only System is MR Safe.
- Respiratory Plate is MR Safe.
- Patient Hand Grips are MR Safe.

FOREHEAD BRACE


NOTE: Device is compatible with Body Pro-Lok™ Bridge and ONEBridge™.

⚠ WARNING

- Ensure device is not too tight for patient comfort.

1. Ensure pins on bridge are in unlatched (*upwards*) position.
2. Align holes and pull accessory towards bridge.
3. Push down pins to latch accessory into place.
4. Attach bridge (*refer to appropriate bridge manual*) with forehead brace (*attached*) at appropriate location.
5. Ensure bridge height adjustment is locked.
6. Lower forehead brace into desired position.

NOTE: Note index marks for desired position.

**SHOULDER BRIDGE 2**

NOTE: Device is compatible with Body Pro-Lok™ Platform and UCT with adapters.


⚠ WARNING

- Ensure device is not too tight for patient comfort.

1. Place shoulder restraint bridge on Body Pro-Lok™ Platform (*refer to appropriate bridge manual*).
2. Adjust height and width of shoulder restraint bridge to approximate height and tighten.
3. Move shoulder restraint bridge in inferior direction until desired shoulder suppression is obtained.

4. Lock shoulder restraint bridge in place using bridge locking knobs.
5. Adjust width and/or height of shoulder restraint bridge as needed.


NOTE: Note all indexed positions of restraints.

RAILS-ONLY SYSTEM

1. Attach Rails-Only System to couchtop at appropriate longitudinal position. Ensure Rails-Only System is secure.

NOTE: • Note all index marks for desired position.
• Note direction of Rails-Only System, with respect to gantry.

RESPIRATORY PLATE


NOTE: Device is compatible with Body Pro-Lok™ Bridge and ONEBridge™.

⚠ WARNING

- Clinical personnel must monitor patient constantly when using respiration suppression accessories to ensure patient can breathe adequately.
- Ensure device is not too tight for patient comfort.

1. Ensure pins on bridge are in unlatched (*upwards*) position.
2. Align holes and pull accessory towards bridge.
3. Push down pins to latch accessory into place.
4. Attach bridge (refer to appropriate bridge manual) with respiratory plate (*attached*) at appropriate location.
5. Lower bridge until device contacts patient and ensure bridge height adjustment is locked.
6. Tighten device until desired respiratory restriction is achieved.

NOTE: Note all indexed positions.


PATIENT HAND GRIPS

NOTE: Device is compatible with Body Pro-Lok™ Platform.

**WARNING**

- Patient must keep hands on handles during device setup, transfer, and motion of platform to prevent injury unless directed otherwise by clinical personnel.

1. Place patient hand grips onto Body Pro-Lok™ Platform at desired position.
2. Turn latching knob to latch handles in place.

REPROCESSING

**WARNING**

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. BODY PRO-LOK, LOK-BAR, ONEBRIDGE AND PRODIGY ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MEDTEC LLC