Body Pro-Lok™ Respiratory Belt

Body Pro-Lok™ System

INSTRUCTIONS FOR USE (IFU)

CE

Towards Gantry

Towards Gantry

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Kızağa Doğru

INDICATIONS FOR USE

The device is indicated to assist in the proper positioning of patients for radiation therapy including electron, photon and proton treatments, and simulation including CT and MR image acquisition.

INTENDED USE

- The device is part of a system intended to immobilize, position and reposition patients undergoing stereotactic radiation therapy including SBRT.
- The device is intended to apply abdominal compression for managing movement and tumor motion during respiration while maintaining maximum comfort to the patient.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Do not use if device appears damaged.
- Inspect device Velcro for damage, wear or contamination that may cause belt to loosen over time.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at www.CQmedical.com.
- Verify patient position with completed setup sheet prior to treatment.
- Ensure patient remains still for duration of setup and treatment.
- Monitor patient constantly when using accessories to ensure patient comfort.
- Ensure clearance is maintained between Body Pro-Lok™ system and all linac or imaging equipment before using system.
- Remove sharp objects from patient and treatment area to avoid puncturing.
- Pressure indication is for reference only.
- Clinical personnel must monitor patient constantly when using respiration suppression accessories to ensure patient can breathe adequately.
- Ensure device is not too tight for patient comfort.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.


MRI SAFETY INFORMATION

⚠ MR Conditional

- Non-clinical testing has demonstrated the Respiratory Belt is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:
 - Static magnetic field of 1.5 T and 3.0 T
 - Maximum spatial field gradient of 3,000 gauss/cm (30 T/m)
 - Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Operating Mode)
- Under the scan conditions defined above, the device is expected to produce a maximum temperature rise of less than 0.5°C after 15 minutes of continuous scanning.
- In non-clinical testing, the image artifact caused by the device extends up to 114mm from the gauge of the Respiratory Belt when imaged with a gradient echo pulse sequence and a 3.0 T MRI System. Image artifact was also observed up to 9mm from the belt clamp (*BPL026 and BPL006 only*).
- Accessory Hand Pump with Pressure Gauge is MR Conditional. Refer to manufacturer’s instructions for conditions of safe use.

1. Attach clamps to side rails.
2. Place Velcro straps loosely through each clamp.
3. Hold respiratory belt firmly in place over desired area.
4. Place Velcro strap over respiratory belt.
5. Place second Velcro strap over respiratory belt, covering first strap if necessary.
6. Attach respiratory gauge and pump to respiratory belt.
7. On gauge, use range A through H for respiratory belt inflation.
8. Inflate respiratory belt until desired respiratory restriction is achieved.
9. Prior to imaging or treatment, remove gauge and pump from respiratory belt.

NOTE: Note all indexed positions and acceptable pressure range on setup sheet.

REPROCESSING

⚠ WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. BODY PRO-LOK AND ONEPLATFORM ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MEDTEC LLC