IPPS™ Overlay

An indexed patient positioning system for use with Philips Ingenia 1.5T and Philips Ingenia 3.0T

INSTRUCTIONS FOR USE (IFU)

IPPS™ Overlay Installation

INSTRUCTION FOR USE (IFU)

INTENDED USE

The device is intended to aid in patient positioning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Do not use if device appears damaged, including but not limited to warpage, cracking, missing components or loose fit.
- Ensure accessories are securely attached in correct position prior to treatment.
- Do not exceed a patient weight of 550 lbs (250 kg).
- Ensure device is secure prior to use.

NOTE:

- If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.
- Overlay weight: 22 lbs (10kg). For total table load rating, including accessories, refer to your system's manual.
- When not in use, store overlay vertical on edge to prevent damage.

MRI SAFETY INFORMATION

MR Safe

- Device is safe for use within MR environment.

ATTACHING OVERLAY TO CRADLE

NOTE: Tools Required (*not provided*): drill and flathead screwdriver.


⚠ WARNING

- Remove cradle from magnet room prior to installing overlay to prevent magnetic interference.


1. Place template (*provided*) in inner channels of cradle, slide to foot end of cradle until template contacts end of channel.

⚠ WARNING

- To obtain proper hole location for overlay attachment and ensure fit of dStream box, template must contact end of channel.


2. Using template, drill three (3) holes a minimum of 1/4" deep into each inside channel of cradle.
3. Remove template, tap center holes using #10-32 tap.
4. Attach locating blocks into inner channels using 10-32 screws.

5. Place overlay on cradle aligning locating blocks with notches on underside of overlay and overlay width retainers (4) into cradle inner channels.


⚠ WARNING

- Grasp overlay at head end and foot end, using two people, to prevent pinched fingers.

STABILIZING LATERAL MOVEMENT OF OVERLAY

1. Loosen locking screws (4) on both sides of overlay.
2. Tighten Adjustment screws (8) until width retainer contacts inside edge of cradle inner channel.
NOTE: Adjust both sides equally so overlay remains centered on cradle.
3. Tightening Locking screws locks width retainer into position.

4. Gently lift overlay off cradle and place back into position to ensure overlay can be installed and removed without excessive effort.

REPROCESSING

⚠ WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

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- NOTE:
- If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.
 - Overlay weight: 22 lbs (10kg). For total table load rating, including accessories, refer to your system's manual.
 - When not in use, store overlay vertical on edge to prevent damage.

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
USING COIL SUPPORTS


1. Position patient as desired.
2. Position coil support bridge in desired location in cradle outside channels.
3. Place anterior coil onto coil support bridge. Using coil support clips, tuck anterior coil underneath clip to secure in place.

**USING HEAD EXTENSION (OPTIONAL)****⚠ WARNING**


- Do not lift overlay using Head Extension, otherwise damage will occur.

NOTE: Usable only if dStream box is removed.

1. Insert attachment pins of Head Extension into head end of overlay. Head Extension should be level with top side of overlay when fully seated.

NOTE: If vertical height adjustment is necessary, remove supports and adjust by adding or removing shims.


2. Position patient as desired.

USING Lok-Bar™

WARNING

- Do not use Lok-Bar™ to secure positioning devices.

1. Snap Lok-Bar™ at appropriate longitudinal position. Ensure Lok-Bar™ is secure.
2. Place positioning accessory onto Lok-Bar™ pins.

WARNING

- Avoid pinch points when overlay is in motion by ensuring patient does not grasp overlay sides.

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1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. IPPS, LOK-BAR, AND PRODIGY ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS.
MEDTEC LLC