
trUpoint ARCH™

Base Lock Replacement Kit

REFERENCE GUIDE


- I. trUpoint ARCH™
- II. Punch (*not included*)
- III. Spring
- IV. Spring Slide
- V. Base Lock
- VI. Flat Blade Screwdriver (*not included*)
- VII. Screw
- VIII. Pin


MRI SAFETY INFORMATION**MR**

MR Safe

- Device is MR safe.
- MR image quality may be compromised if the area of interest is relatively close to the device. Optimize MR imaging parameters prior to treating patient. Image artifact extends approximately 12mm from device when scanned in nonclinical testing using a gradient echo sequence in a 3.0T Siemens Skyra MR system with built-in body coil.

REMOVING BASE LOCK ASSEMBLY


1. Remove screw with flat blade screwdriver.
2. Push pin out using punch.
3. Remove base lock, spring slide and spring and discard.

ATTACHING BASE LOCK ASSEMBLY

1. Place spring, spring slide and base lock into trUpoint ARCH™ and hold in place.
2. Push pin firmly into trUpoint ARCH™ using punch.

3. Insert and tighten screw using flat blade screwdriver.

-  **WARNING**
- *Ensure screw is engaged fully to block pin hole.*
 - *Ensure base locks are secure and trUpoint ARCH™ is securely attached to baseplate prior to setup and treatment of patient.*
 - *Inspect device periodically for signs of damage and general wear.*
 - *Do not use if device appears damaged.*
 - *When not in use, store on trUpoint ARCH™ Storage Plate to prevent damage (refer to trUpoint ARCH™ Storage Plate instruction).*

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. TRUPOINT ARCH IS A TRADEMARK OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MEDTEC LLC