ONEBridge™

Body Pro-Lok™ System

INSTRUCTIONS FOR USE (IFU)

CE

 Towards Gantry	Towards Gantry باتجاه القطار 朝向扫描架 Prema kućištu Směrem k nosné konstrukci Mod ophængning Naar rijbrug Portaali suunas Pöytää kohti Couvercle de portique Zum Gestell hin Προς πλαίσιο A gantry irányába Verso il gantry 構台方向 갠트리 방향 Gentrija virzienā Link rémo Mot gantry W kierunku suwnicy Em direção à estrutura Em direção à torre Spre suportul de scanare К гентрик гентри Prema postolju Smerom k portálu Proti podstavku Hacia la pasarela Mot gantry ไปทางโครงสร้างลิ้นชัก Kızağa Doğru
--	---

INTENDED USE

The device is intended to aid in support and positioning of patients during MR, radiological and other procedures.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Do not use if device appears damaged.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at www.CQmedical.com.
- Ensure patient is set up in correct position before beginning treatment.
- Ensure patient remains still for duration of setup and treatment.
- Monitor patient constantly when using accessories to ensure patient comfort.
- Ensure clearance is maintained between Body Pro-Lok™ system and all linac or imaging equipment before using system.
- Do not attempt to move or carry platform by accessory bridges.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.


MRI SAFETY INFORMATION

- Device is MR safe.


MR Safe

NOTE: MR image quality may be compromised if the area of interest is relatively close to the bridge base and bridge side. Optimize MR imaging parameters prior to treating patient. Image artifact extends approximately 11mm from bridge docking base and bridge sleeve and 15mm from bridge arm when scanned in nonclinical testing using a gradient echo sequence in a 3.0T Siemens Skyra MR system with built-in body coil.

1. Loosen lateral adjustment levers.
2. Raise bridge to highest position and lock height adjustment levers.
3. Ensure docking levers are in unlocked position.
4. Place bridge onto patient support at appropriate location.
5. Lock bridge docking lever, ensuring bridge is securely affixed.
6. Attach accessory, if applicable.
7. Adjust bridge height and lock height adjustment lever at desired setting.
8. Adjust bridge laterally, as needed.
9. Tighten lateral adjustment levers.

**QUICK RELEASE**

1. Unlock docking levers on one side.
2. Lift released side of bridge clear of patient.

REPROCESSING

WARNING

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*

1. If necessary, clean surface by removing visual contamination with 70% alcohol or Theracide (PDI Sani-Cloth Plus) wipes. If visual contamination cannot be removed, repeat cleaning steps and if necessary, discard device.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

STORAGE

NOTE: Store bridge with lateral adjustment levers and docking levers in unlocked position.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. BODY PRO-LOK, ONEBRIDGE AND ONEPLATFORM ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MEDTEC LLC