Universal Couchtop Accessories

INSTRUCTIONS FOR USE (IFU)

CQ MEDICALTM

CE

INTENDED USE

The device is indicated to aid in supporting and positioning adult and pediatric patients undergoing radiation therapy including electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Do not use if device appears damaged.
- Verify all angles of treatment, attenuation characteristics and WET values prior to treating patients.


- NOTE:
- If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.
 - Devices are compatible with Universal Couchtop™ and ONEPlatform™.

BRIDGE ADAPTER

NOTE: Body Pro-Lok™ bridges with an adapter can be used with Universal Couchtop™ and ONEPlatform™.

1. Attach bridge adapter to clamps.


NOTE: When using adapter on Universal Couchtop™ or ONEPlatform™, bridge must be set one index mark higher (5mm) than on the Body Pro-Lok™ platform.

RESPIRATORY BELT ADAPTER

NOTE: For use with Respiratory Belt.

1. Attach respiratory belt adapter to patient support.
2. Attach clamp to respiratory belt adapter.


ACCESSORY RAIL ADAPTER

⚠ WARNING

- Do not exceed load of 44 lbs (20 kg).
- Ensure all accessories are securely attached.

1. Place accessory rail adapter onto patient support at desired position.
2. Turn locking knob to secure rail in place.

PATIENT HAND GRIPS

⚠ WARNING

- Patient must keep hands on hand grips during device setup to prevent injury unless directed otherwise by clinical personnel.


MRI SAFETY INFORMATION

- MTIL662501, MTIL662502: Device is MR safe.

MR Safe

1. Place patient hand grips onto patient support at desired position.
2. Turn locking knob to secure handles in place.

REPROCESSING

⚠ WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol.
- NOTE: Color may transfer to wipe during normal cleaning.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2023 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. BODY PRO-LOK, ONEPLATFORM AND UNIVERSAL COUCHTOP ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
**MEDTEC LLC**