# Bite Cup

*For use with trUpoint ARCH™*

INSTRUCTIONS FOR USE (IFU)

---

**INTENDED USE**

Device is intended to be used for immobilization, positioning and repositioning during Stereotactic Radiotherapy (SRT) and Stereotactic Radiosurgery (SRS) in all parts of the brain, head, and neck during external beam radiation therapy. Device is also used to immobilize and position the head during image acquisition to support treatment planning including in Computed Tomography (CT) and Magnetic Resonance (MR) imaging systems.

**CAUTION**

Federal (United States) law restricts this device to sale by or on the order of a physician.

**⚠ WARNING**

- Do not use if device appears damaged.
- Inspect single patient components for damage or contamination prior to use. If evident, dispose of as infectious waste and create new components.
- trUpoint ARCH™ is for use with CQ Medical specified accessories only.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at [www.CQmedical.com](http://www.CQmedical.com).
- Verify patient position with verification image and completed setup sheet prior to treatment.
- Prior to patient use, ensure all adjustments on head frame are secure and baseplate is stable on surface.
- A confirmation image is recommended to ensure patient and target are in the desired location before treatment is delivered.
- Ensure clearance is maintained between trUpoint ARCH™ and all linac or imaging equipment.

**NOTE:** If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

**MRI SAFETY INFORMATION**

- Device is MR safe.

MR Safe

**NOTE:** MR image quality may be compromised if the area of interest is relatively close to the bite cup and bite tray filled with putty, nasion cushion and arch. Optimize MR imaging parameters prior to treating patient. Image artifact extends approximately 19mm from bite cup and bite tray filled with putty, 16mm from nasion cushion and 12mm from the arch when scanned in nonclinical testing using a gradient echo sequence in a 3.0T Siemens Skyra MR system with built-in body coil.


**USING BITE ASSEMBLY**

**⚠ WARNING**

- Bite cup is intended for single patient throughout setup and treatment cycle.
- Ensure bite tray is fully secure in bite cup prior to treatment.

**INITIAL BITE CUP SET-UP**


1. Unlock black lever on bite assembly. Loosen small knob on bite assembly and medium knob on bite assembly arm prior to initial set-up.


2. With bite tray in patient's mouth, hold bite cup directly over bite tray stem. Slide bite assembly arm toward bite cup to obtain desired linear position and adjust bite assembly arm to desired angle. Tighten small knob on bite assembly and medium knob on assembly arm to secure.
3. Adjust bite cup to desired vertical position and lock black lever on bite assembly to secure.

NOTE:


- Position bite cup as close to patient's lips as possible without touching.
- Ensure arch of bite tray stem is visible in the top two holes of bite cup.


4. Insert putty (refer to putty instructions for additional warnings) into bite cup around bite tray stem. Ensure putty fills the bite cup completely, indicated by putty protruding from holes in bite cup. Allow putty to fully harden.

#### ADJUSTING BITE ASSEMBLY

**ADJUSTING BITE CUP****DISPOSAL****⚠ WARNING**

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*
- *Dispose of single patient components as infectious waste.*

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. TRUPOINT ARCH IS A TRADEMARK OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
**MEDTEC LLC**