Nasion Cushion

For use with trUpoint ARCH™

INSTRUCTIONS FOR USE (IFU)

CE

INTENDED USE

Device is intended to be used for immobilization, positioning and repositioning during Stereotactic Radiotherapy (SRT) and Stereotactic Radiosurgery (SRS) in all parts of the brain, head, and neck during external beam radiation therapy. Device is also used to immobilize and position the head during image acquisition to support treatment planning including in Computed Tomography (CT) and Magnetic Resonance (MR) imaging systems.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Do not use if device appears damaged.
- Inspect single patient components for damage or contamination prior to use. If evident, dispose of as infectious waste and create new components.
- trUpoint ARCH™ is for use with CQ Medical specified accessories only.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at www.CQmedical.com.
- Verify patient position with verification image and completed setup sheet prior to treatment.
- Prior to patient use, ensure all adjustments on head frame are secure and baseplate is stable on surface.
- A confirmation image is recommended to ensure patient and target are in the desired location before treatment is delivered.
- Ensure clearance is maintained between trUpoint ARCH™ and all linac or imaging equipment.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

MRI SAFETY INFORMATION

- Device is MR safe.


MR Safe

NOTE: MR image quality may be compromised if the area of interest is relatively close to the bite cup and bite tray filled with putty, nasion cushion and arch. Optimize MR imaging parameters prior to treating patient. Image artifact extends approximately 19mm from bite cup and bite tray filled with putty, 16mm from nasion cushion and 12mm from the arch when scanned in nonclinical testing using a gradient echo sequence in a 3.0T Siemens Skyra MR system with built-in body coil.

USING NASION ASSEMBLY**⚠ WARNING**


- Nasion cushion is intended for single patient throughout setup and treatment cycle.

1. Attach nasion cushion to nasion stem.
2. Position nasion assembly in highest position and tighten small knob.

3. Loosen small knob on nasion assembly.
4. Position nasion cushion at lower arch of patient's nose bridge junction, making vertical adjustment last. Tighten small knob to secure.

ADJUSTING NASION ASSEMBLY

DISPOSAL

⚠ WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.
- Dispose of single patient components as infectious waste.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. TRUPOINT ARCH IS A TRADEMARK OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MEDTEC LLC