# Monarch™ Overhead Arm Positioner

*For use with C-Qual™*

INSTRUCTIONS FOR USE (IFU)

---

CE

**INTENDED USE**

The intended use of this device is to support and aid in positioning a patient during radiologic and other medical procedures.

**CAUTION**

Federal (United States) law restricts this device to sale by or on the order of a physician.



**⚠ WARNING**

- Do not use if device appears damaged.
- Do not allow patient to reposition themselves.
- Do not allow patient to reposition using hand grip as these are not intended to be load bearing.
- Do not reposition device with patient on it.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at [www.CQmedical.com](http://www.CQmedical.com).
- Verify patient position with completed setup sheet prior to treatment.
- Do not exceed 200 lbs (91 kg) on wings.
- Ensure device is secure before positioning patient.
- Ensure device/patient does not move for duration of set-up and treatment.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Verification image should be taken to confirm patient position prior to treatment.
- Ensure device is secure prior to use.
- Do not use device with unapproved accessories or devices.


**NOTE:** If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

 MR Conditional	<ul style="list-style-type: none"> <li>• Non-clinical testing and scientific rationale has demonstrated the C-Qual M™ Breastboard with Monarch™ Overhead Arm Positioner is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:                     <ul style="list-style-type: none"> <li>• Static magnetic field of 1.5 T and 3.0 T</li> <li>• Maximum spatial field gradient of 500 gauss/cm (5 T/m)</li> <li>• Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Operating Mode)</li> </ul> </li> <li>• Under the scan conditions defined above, the device is expected to produce a maximum temperature rise of less than 1° C after 15 minutes of continuous scanning.</li> <li>• In non-clinical testing, image artifacts are present outside a typical area of treatment. Near the superior locking knob, the artifact extends up to 124mm when imaged with a gradient echo pulse sequence and a 3.0 T MRI system. Image artifacts were also found at other locations, but with less severity.</li> </ul>
---	--


1. Remove arm support and detachable hardware.

	<p>I. Arm Supports and Detachable Hardware</p>
	

2. Install Monarch™ brackets with provided hardware.


	<p>I. Brackets                  II. Bracket Hardware (M6)                  III. Bracket Hardware (M5)</p>
---	---

3. Install crossbar with provided hardware.


NOTE: Ensure device is properly assembled and secure.

4. Attach guide bar assembly by tightening thumbscrew.


5. Insert locking knob assemblies as shown.

6. Turn ellipse knobs 90 degrees to secure.


NOTE: Ensure upper board is in locked position.


7. Slide Monarch™ onto C-Qual™ breastboard.
8. Secure locking knobs.
9. Place desired head support into recessed area of Monarch™.


NOTE: Ensure to not overtighten locking knobs.

10. Adjust C-Qual™ as needed (refer to C-Qual™ instruction).
11. Adjust hand grip as needed.


Adjust Monarch™ Individual Hand Grips as needed.

12. Adjust Monarch™ as needed.

**REMOVING MONARCH™ OVERHEAD ARM POSITIONER FROM C-QUAL™**

NOTE: • Prior to removing Monarch™, lock upper board to lower board.  
 • Monarch™ not shown.

1. Unlock locking knobs on Monarch™.
2. Slide Monarch™ clear of locking brackets and lift off of breastboard.

**REPROCESSING**

**⚠ WARNING**

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.
- Inspect device for damage such as cracks or breakage. If damage is evident, discontinue use and contact CQ Medical for further instructions as needed.

1. Clean device with a mild, non-abrasive cleaner as needed.

**MAINTENANCE**

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. C-QUAL AND MONARCH ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS.
MEDTEC LLC