Body Pro-Lok™ Respiratory Plate with ComfortCare™ Cushion

Body Pro-Lok™ System

INSTRUCTIONS FOR USE (IFU)

CE

INDICATIONS FOR USE

The device is indicated to assist in the proper positioning of patients for radiation therapy including electron, photon and proton treatments, and simulation including CT and MR image acquisition.

INTENDED USE

- The device is part of a system intended to immobilize, position and reposition patients undergoing stereotactic radiation therapy including SBRT.
- The device is intended to apply abdominal compression for managing movement and tumor motion during respiration while maintaining maximum comfort to the patient.

MRI SAFETY INFORMATION

MR Safe

- Respiratory Plate is MR Safe.

MR Conditional

- Non-clinical testing and scientific rationale has demonstrated that the Respiratory Cushion is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:
 - Static magnetic field of 1.5 T and 3.0 T
 - Maximum spatial field gradient of 3,000 gauss/cm (30 T/m)
 - Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Operating Mode)
- Under the scan conditions defined above, the device is expected to produce a maximum temperature rise of less than 1° C after 15 minutes of continuous scanning.
- In non-clinical testing, no image artifact was caused by the device when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

MR Conditional

- Accessory Hand Pump with Pressure Gauge is MR Conditional. Refer to manufacturer's instructions for conditions of safe use.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.


WARNING

- Do not use if device appears damaged.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at www.CQmedical.com.
- Verify patient position with completed setup sheet prior to treatment.
- Monitor patient constantly when using accessories to ensure patient comfort.
- Ensure clearance is maintained between Body Pro-Lok™ system and all linac or imaging equipment before using system.
- Remove sharp objects from patient and treatment area to avoid puncturing.
- Pressure indication is for reference only.
- Clinical personnel must monitor patient constantly when using respiration suppression accessories to ensure patient can breathe adequately.
- Ensure device is not too tight for patient comfort.
- Inflate ComfortCare™ cushion using CQ Medical supplied hand pump only.
- Ensure patient remains still for duration of setup and treatment.

- NOTE:**
- If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.
 - Device is compatible with Body Pro-Lok™ Bridge and ONEBridge™.

1. Ensure cushion is deflated prior to setting up patient.
2. Attach ComfortCare™ cushion to respiratory plate by placing the valve through the round hole and edges of plate into pockets on cushion.
3. Ensure pins on bridge are in unlatched (upwards) position.
4. Align holes and pull respiratory plate towards bridge.
5. Push down pins to latch respiratory plate into place.
6. Attach bridge (refer to appropriate bridge manual) with respiratory plate (attached) at appropriate location.
7. Lower bridge until device contacts patient and ensure bridge height adjustment is locked.
8. Lower plate until desired respiratory restriction is achieved.
9. Inflate ComfortCare™ cushion to distribute pressure and improve patient comfort.
10. Prior to imaging or treatment, remove pump and gauge from cushion.

- NOTE:**
- Use range 1 through 8 on respiratory gauge for ComfortCare™ cushion inflation.
 - Record all patient specific settings on setup sheet.

REPROCESSING **WARNING**

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*

1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol.

MAINTENANCE

- NOTE:
- Inspect device prior to use for signs of damage and general wear.
 - ComfortCare™ cushion is designed for a maximum of 200 uses.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. BODY PRO-LOK, COMFORTCARE AND ONEBRIDGE ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MEDTEC LLC