SnapFrame™ Type-S™

REFERENCE GUIDE

CIVCO

CE

INTENDED USE

The device is indicated to aid in supporting, positioning, and/or immobilization of adult and pediatric patients undergoing radiation therapy of the head, brain, neck, and spine including radiosurgery and electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

MRI SAFETY INFORMATION

- Devices are MR Safe.

MR Safe

⚠ WARNING

- Do not use if device appears damaged.
- Device should be used for a single patient throughout setup and treatment cycle.
- Ensure baseplate is secure on treatment platform.
- Cutting thermoplastic masks may reduce the rigidity of the device. If cutting is necessary, ensure rigidity of mask is maintained by adding thermoplastic strips. Thermoplastic strips are not included with masks but can be purchased by contacting your CIVCO representative.
- Ensure thermoplastic is properly aligned over patient before forming.
- Ensure nose mouth portal is aligned during molding to allow patient to breathe from nose and mouth during future use.
- Ensure thermoplastic is formed properly around treatment area.
- Form thermoplastic to avoid contact with tracheotomy tube.
- Ensure thermoplastic is properly attached using clips and pins.
- Verify all angles of treatment, attenuation characteristics and WET values prior to treating patients. Refer to Technical Data Sheet at WWW.CIVCORT.COM for WET values.
- Thermoplastic material may cause skin irritation in individuals with sensitive skin.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

USING THERMOPLASTICS

1. Setup waterbath with 2–3" of water and heat (refer to table for appropriate temperature).

Thermoplastic Type	Waterbath Temperature	Heating Time
FirmFit™	70°C (158°F)	1 min

2. Place thermoplastic into bath for directed time.
- NOTE: Improper heating time may cause difficult application.


3. Remove thermoplastic from waterbath.
4. Quickly remove excess water from thermoplastic (less than five seconds).
5. Position thermoplastic over treatment area.

- NOTE:
- Overstretching of thermoplastic material may affect rigidity of mask. After pulling thermoplastic mask, ensure rigidity of mask is maintained.
 - Some thermoplastic masks may require assistance in pulling of mask.


6. Gently pull sides of frame outward and down to attachment base. Secure pins.

NOTE: Open Face Thermoplastic Only:

- Ensure thermoplastic is secure at forehead and chin locations when pulling to attachment base.
- Ensure desired anatomical locations are free of thermoplastic material as required by camera system region of interest.


7. With sides of frame attached to base, pull top section of frame outward and down to attachment base, while securing thermoplastic at chin. Secure pin and frame clips.

8. Head, Neck and Shoulder Only: Pull shoulder frame outward and down to attachment base. Secure pins.

NOTE: Head, Neck and Shoulder Thermoplastic Only:

- Ensure thermoplastic is secure at chest area when pulling to attachment base.

9. Gently mold material to contours of treatment area.
10. Allow thermoplastic to cool completely before releasing from attachment base.

REPROCESSING

WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. If necessary, clean surface of thermoplastic and/or frames and pins by removing gross visual contaminants with alcohol or warm soapy water. If gross visual contamination cannot be removed repeat cleaning steps and if necessary, discard device.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2021 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. FIRMFIT, SNAPFRAME AND TYPE-S ARE TRADEMARKS OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.