ProForm™ MR Overlay

INSTRUCTIONS FOR USE (IFU)

INTENDED USE

The device is indicated to aid in supporting, positioning, and/or immobilization of adult and pediatric patients undergoing radiation therapy of the head, brain, neck, and spine including radiosurgery and electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Do not use if device appears damaged.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Prior to treatment, ensure device and accessories are secure.
- Lok-Bar™ is a positioning device only, do not use to secure patient.
- Do not use Lok-Bar™ to secure positioning devices.
- Ensure device/patient does not move for duration of set-up and treatment.
- Do not reposition device with patient on it.
- Do not allow patient to reposition themselves, such as arching back.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

MRI SAFETY INFORMATION


MR Safe

- Device is MR safe.

USING PROFORM™ MR OVERLAY

NOTE: Ensure overlay is fully supported by couchtop.

1. Place overlay onto Lok-Bars (2) (refer to Lok-Bar™ instructions).

2. Position patient as dictated by setup protocol.
3. Referring to AccuForm™ Cushion and ProForm™ Thermoplastic instructions, form AccuForm™ Cushion and ProForm™ Thermoplastic as needed.

REPROCESSING**⚠ WARNING**

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Using a lint-free cloth wetted with tap water, wipe away gross soil.
2. Using PDI Sani-Cloth Plus wipes, wipe articles again to further remove soil.
 - a. Repeat this step one more time for a total of two wiping steps.
3. Using a lint-free cloth wetted with purified water, further wipe away residuals.
 - a. Repeat this step one more time for a total of two wiping steps.
4. Allow articles to air dry.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. LOK-BAR AND PROFORM ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.