Solstat™ Immobilization System

INSTRUCTIONS FOR USE (IFU)

CE

INDICATIONS FOR USE / INTENDED USE

The device is indicated to position and/or immobilize adult and pediatric patients undergoing radiation therapy of the head, brain, and neck, including Surface Guided Radiation Therapy (SGRT) and electron, photon, and proton treatments. The device is also used during image acquisition, including Computed Tomography (CT) Magnetic Resonance (MR) Imaging, to support treatment planning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Verify all angles of treatment, attenuation characteristics and WET values prior to treating patients.
- Product shape may have slight variations; the same device should be used for a patient throughout setup and treatment cycles.
- Do not use if device appears damaged.
- Ensure accessory device is properly located.
- Do not use without accessory cushion or headrest.
- Device is for use with CQ Medical specified accessories only.
- Ensure device is secure prior to use.
- Verification image should be taken to confirm patient position prior to treatment.
- Do not use excessive force on device. Injury or device damage may occur.
- Ensure device/patient does not move for duration of set-up and treatment.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.


MRI SAFETY INFORMATION

MR Safe


- Device is MR safe.

USING SOLSTAT™ IMMOBILIZATION SYSTEM

NOTE: Ensure pins are unlatched.

1. Attach head fixation to extension.
2. Use all available pins to secure head fixation to extension.

NOTE: If using AccuForm™ cushion, proceed to steps 3-9. If using headrest, refer to manual provided with headrest.

3. Prepare AccuForm™ cushion - press spray bottle nozzle slightly into cushion and squeeze a shot of water into cushion at one inch intervals across surface to ensure even dampening. Dab cushion with dry towel to remove excess water.
4. Flatten AccuForm™ cushion to a uniform thickness.
5. Hold cushion over device opening. Align superior edge of cushion with superior edge of bowl.

⚠ WARNING

- Ensure cushion is positioned below edge of bowl to allow attachment of thermoplastic.

6. Lower patient's head into head support.
7. Allow cushion to set (*harden*) approximately 20 minutes.

NOTE: Instruct patient not to move while cushion is setting.

8. Remove patient from cushion without putting pressure on cushion to avoid deforming cushion.
9. Allow cushion to air dry completely (*approximately 24 hours*).

⚠ WARNING

- Do not use heat to dry cushion.

10. Attach thermoplastic (*refer to thermoplastic manual*).

REPROCESSING**WARNING**

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*

1. Clean surface of device by removing visual contaminants with common germicidal or antiseptic wipe, such as quaternary ammonium or isopropyl alcohol wipes. If visual contamination cannot be removed, repeat cleaning steps and if necessary, discontinue use.

MAINTENANCE

NOTE: Inspect for damage such as cracks, breakage, or warpage.

STORAGE**WARNING**

- *Store device on a flat surface. Do not stack.*

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. SOLSTAT AND SOLSTICE ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS.