Headrest

*For use with Solstat™ Immobilization System and Solstice™ SRS
Immobilization System*

INSTRUCTIONS FOR USE (IFU)

INDICATIONS FOR USE / INTENDED USE

The intended use of this device is to immobilize the head during radiation therapy and to support and position during an MRI.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

MRI SAFETY INFORMATION

MR Safe

- Device is MR safe.

WARNING


- *Verify all angles of treatment, attenuation characteristics and WET values prior to treating patients.*
- *Product shape may have slight variations; the same device should be used for a patient throughout setup and treatment cycles.*
- *Do not use if device appears damaged.*
- *Ensure accessory device is properly located.*
- *Ensure device is secure prior to use.*
- *Verification image should be taken to confirm patient position prior to treatment.*
- *Ensure device/patient does not move for duration of set-up and treatment.*
- *Do not use excessive force on device. Injury or device damage may occur.*
- *Device is for use with CQ Medical specified accessories only.*

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

USING HEADRESTS

1. Place headrest onto device.

NOTE: Ensure notches are placed into locating features.

2. Position patient on headrest.

REPROCESSING

WARNING

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*

1. Clean surface of device by removing visual contaminants with common germicidal or antiseptic wipe, such as alcohol. If visual contamination cannot be removed, repeat cleaning steps and if necessary, discontinue use.

MAINTENANCE

NOTE: Inspect for damage such as cracks, breakage, or warpage.

STORAGE

WARNING

- *Store headrests on a flat surface. Do not stack.*

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. SOLSTAT AND SOLSTICE ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
**MEDTEC LLC**