Combifix™

INSTRUCTIONS FOR USE (IFU)

INTENDED USE

The device is indicated to aid in supporting and positioning adult and pediatric patients undergoing radiation therapy including electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- *Verify all angles of treatment, attenuation characteristics and WET values prior to treating patients.*
- *Do not use if device appears damaged.*
- *Ensure all accessories and parts are securely placed prior to treatment.*
- *Do not reposition device with patient on it.*
- *When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at www.CQmedical.com.*
- *Verify patient position with completed setup sheet prior to treatment.*
- *Device is for use with CQ Medical specified accessories only.*
- *Ensure device is secure prior to use.*
- *Do not use excessive force on device. Injury or device damage may occur.*
- *Use MR compatible accessories only.*
- *Ensure device/patient does not move for duration of set-up and treatment.*

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.


MRI SAFETY INFORMATION

- MR Conditional
- Non-clinical testing and scientific rationale has demonstrated that the Combifix™ is MR Conditional. A patient with this device can be safely scanned in a MR System meeting the following conditions:
 - Static magnetic field of 1.5 T and 3.0 T
 - Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Operating Mode)
 - Under the scan conditions defined above, the device is expected to produce a maximum temperature rise of less than 0.6°C after 15 minutes of continuous scanning.


DIRECTIONS FOR USE

1. Place Combifix™ board on Lok-Bar™ pins (See Lok-Bar™ instructions for Lok-Bar™ use).


I. Optional


2. Place Feetfix™ cushion holder on Combifix™ board, ensuring that Feetfix™ knobs fit into slots on Combifix™ board.

3. Place Feetfix™ cushion in Feetfix™ holder and Kneefix™ in appropriate location with desired elevation blocks (optional).


4. Assist patient in lying down on table and placing legs over Kneefix™ and placing feet in Feetfix™.
5. Adjust Combifix™ as needed for comfort of patient and desired treatment conditions.

6. Record all positions for setup as needed on setup sheet.

NOTE: If desired, the Feetfix™ holder can be used directly on a table top with a 2-pin Lok-Bar™. To do so, remove the knobs by removing the screws as shown.

REPROCESSING

WARNING

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*

NOTE:

- If device becomes contaminated or can no longer be cleaned, discontinue use and replace.
- Device must be inspected before and after use for damage and wear.

1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. COMBIFIX, FEETFIX, KNEEFIX AND LOK-BAR ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MEDTEC LLC