Clam-Lok™ Cushion

INSTRUCTIONS FOR USE (IFU)

CE

INTENDED USE

Device is intended to immobilize a patient during radiation therapy and is used to aid in the support and positioning of patients during an MRI.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- *Verify all angles of treatment and attenuation characteristics prior to treating patients.*
- *Do not use if device appears damaged.*
- *Ensure Clam-Lok™ cushion is deflated to proper vacuum prior to each use.*
- *Ensure device is properly attached using bridge latching mechanisms.*
- *When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at www.CQmedical.com.*
- *Verify patient position with completed setup sheet prior to treatment.*
- *Ensure patient remains still for duration of setup and treatment.*
- *Remove sharp objects from patient and treatment area to avoid puncturing.*
- *Ensure Clam-Lok™ cushion is not over-inflated to avoid bursting or difficulty in positioning patient.*
- *Material may cause skin irritation in individuals with sensitive skin.*

- NOTE:**
- If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.
 - Monitor patient constantly when using accessories to ensure patient comfort.

MRI SAFETY INFORMATION

MR
Conditional


- Clam-Lok™ Cushion is MR conditional. Nonclinical testing has demonstrated that the Clam-Lok™ Cushion can be scanned safely at 1.5T and 3.0T environments. In nonclinical testing, the Clam-Lok™ Cushion valve stem created marginal artifacts at 1.5T and 3.0T environments.
 - Static magnetic field of 1.5 T and 3.0 T
 - Maximum spatial field gradient of 2,000 gauss/cm (20 T/m)
 - Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Operating Mode)
- Under the scan conditions defined above, the device is expected to produce a maximum temperature rise of less than 0.4°C after 15 minutes of continuous scanning.
- In non-clinical testing, the image artifact caused by the device extends up to 82mm from the Clam-Lok™ valve with stem when imaged with a spin echo pulse sequence and 3.0 T MRI System.

CLAM-LOK™ CUSHION

- NOTE:**
- Device is compatible with Body Pro-Lok™ Bridge and ONEBridge™.
 - Ensure integrity of Clam-Lok™ cushion is not compromised during daily handling.
 - Clam-Lok™ cushion should be inflated prior to use.

1. Ensure pins on bridge are in unlatched (*upwards*) position.
2. Align holes and pull accessory towards bridge.
3. Push down pins to latch accessory into place with bridge at maximum height.
4. Attach bridge (*refer to appropriate bridge manual*) with Clam-Lok™ Cushion (*attached*) at appropriate location.
5. Adjust bridge height to desired setting.
6. Ensure bridge height adjustment is locked.

7. Inflate Clam-Lok™ cushion to moldable consistency using Vac-Lok™ Vacuum Compressor (refer to Vac-Lok™ Vacuum Compressor instruction) or place female valve connector onto cushion valve to let air into Clam-Lok™ cushion without using Vacuum Compressor.
8. Form Clam-Lok™ cushion around appropriate location and evacuate air from Clam-Lok™ cushions. (See Vac-Lok™ instruction)
9. Ensure device is not too tight for patient comfort.

REPROCESSING

WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. If necessary, clean surface by removing visual contamination with common germicidal wipe, such as alcohol. If visual contamination cannot be removed repeat cleaning steps and if necessary, discard device.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. BODY PRO-LOK, CLAM-LOK, ONEBRIDGE AND VAC-LOK ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
**MEDTEC LLC**