Monarch™ Individual Hand Grips

For use with C-Qual™

REFERENCE GUIDE

CIVCO
Radiotherapy™

CE

INTENDED USE

- The device is intended to support and aid in positioning a patient during radiologic and other medical procedures.
- Devices are used to aid in the support and positioning of patients during an MRI.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Do not use if device appears damaged.
- Do not use excessive force on device. Injury or device damage may occur.
- Do not reposition device with patient on it.
- Ensure device is secure before positioning patient.
- Ensure device/patient does not move for duration of set-up and treatment.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at WWW.CIVCORT.COM.
- Verify patient position with completed setup sheet prior to treatment.
- Do not use device with unapproved accessories or devices.
- Do not allow patient to reposition using hand grip as these are not intended to be load bearing.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.


MRI SAFETY INFORMATION


- Non-clinical testing has demonstrated the Monarch™ Individual Hand Grips are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:
 - Static magnetic field of 1.5 T and 3.0 T
 - Maximum spatial field gradient of 500 gauss/cm (5 T/m)
 - Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Operating Mode)
- Under the scan conditions defined above, the device is expected to produce a maximum temperature rise of less than 1° C after 15 minutes of continuous scanning.

INSTALLING MONARCH™ INDIVIDUAL HAND GRIPS ON MONARCH™ PATIENT SUPPORT

1. If necessary, remove existing Hand Grip Tracks from Monarch™ Patient Support and install new Hand Grip Tracks with provided hardware.

2. Adjust Monarch™ Individual Hand Grips as needed.

REPROCESSING

WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

- If necessary, clean surface by removing visual contamination with common germicidal wipe, such as alcohol. If visual contamination cannot be removed repeat cleaning steps and if necessary, discard device.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2022 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. CIVCO RADIOTHERAPY AND CIVCO RT ARE TRADEMARKS OF MEDTEC LLC. MONARCH IS A TRADEMARK OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.