# Solstice™ MR Immobilization System

---

REFERENCE GUIDE

**CIVCO**

**INDICATIONS FOR USE / INTENDED USE**

The device is indicated to position and/or immobilize adult and pediatric patients undergoing radiation therapy of the head, brain, and neck, including Surface Guided Radiation Therapy (SGRT) and electron, photon, and proton treatments. The device is also used during image acquisition, including Computed Tomography (CT) Magnetic Resonance (MR) Imaging, to support treatment planning.

**CAUTION**

Federal (United States) law restricts this device to sale by or on the order of a physician.

**⚠ WARNING**

- Verify all angles of treatment, attenuation characteristics and WET values prior to treating patients.
- Product shape may have slight variations; the same device must be used for a patient throughout setup and treatment cycles.
- Do not use if device appears damaged.
- Ensure accessory device is properly located.
- Device is for use with CIVCO specified accessories only.
- Ensure device is secure prior to use.
- Verification image should be taken to confirm patient position prior to treatment.
- Do not use excessive force on device. Injury or device damage may occur.
- Ensure device/patient does not move for duration of set-up and treatment.
- Do not use without accessory cushion or headrest.

**MRI SAFETY INFORMATION**


MR Safe


- Device is MR safe.

**USING SOLSTICE™ MR IMMOBILIZATION SYSTEM**

NOTE: Ensure pins are unlatched.


1. Attach head fixation to extension.
2. Use all available pins to secure head fixation to extension.

**⚠ WARNING**

- Bowl indexing does not relate to patient isocenter.

3. For initial setup, position indexing location at 5.

NOTE: If using AccuForm™ cushion, proceed to steps 4-10. If using headrest, refer to manual provided with headrest.

4. Prepare AccuForm™ cushion - press spray bottle nozzle slightly into cushion and squeeze a shot of water into cushion at one inch intervals across surface to ensure even dampening. Dab cushion with dry towel to remove excess water.
5. Flatten AccuForm™ cushion to a uniform thickness.
6. Hold cushion over device opening. Align superior edge of cushion with superior edge of bowl.
7. Lower patient's head into head support.

 **WARNING**

- Position cushion below thermoplastic attachment edge of bowl, ensuring inferior edge of cushion at neck area does not limit pitch functionality.

8. Allow cushion to set (*harden*) approximately 20 minutes.  
NOTE: Instruct patient not to move while cushion is setting.
9. Remove patient from cushion without putting pressure on cushion to avoid deforming cushion.
10. Allow cushion to air dry completely (*approximately 24 hours*).

 **WARNING**

- Do not use heat to dry cushion.

11. Attach thermoplastic (*refer to thermoplastic manual*).
12. Adjust bowl indexing location as needed.

- NOTE:
- Rotate screw clockwise to lower pitch.
  - Rotate screw counter-clockwise to raise pitch.

**REPROCESSING** **WARNING**

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Clean surface of device by removing visual contaminants with common germicidal or antiseptic wipe, such as quaternary ammonium or isopropyl alcohol wipes. If visual contamination cannot be removed, repeat cleaning steps and if necessary, discontinue use.

**MAINTENANCE**

- NOTE:
- Inspect for damage such as cracks, breakage, or warpage.
  - User should not service device. If service is needed, contact CIVCO.

**STORAGE** **WARNING**

- Do not stack device.
- Store device on flat surface.

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COPYRIGHT © 2020 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. SOLSTICE IS A TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.