# IPPS™ Overlay

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**INTENDED USE**

Devices are used to aid in the support and positioning of patients during an MRI.

**CAUTION**

Federal (United States) law restricts this device to sale by or on the order of a physician.

**⚠ WARNING**

- Do not use if device appears damaged.
- Ensure all accessories are securely attached in correct position prior to imaging.
- Do not apply a load exceeding 550 lbs (250 kg).
- Ensure device is secure prior to use.
- Store overlay in upright position when not in use to prevent damage to overlay supports.
- Use MR compatible accessories only.
- Avoid scanning in first level mode.
- Avoid scanning sedated patients.
- Ensure patient is informed to use squeeze ball in case of heating sensation.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

**MRI SAFETY INFORMATION**

- Device is MR safe.


MR Safe

**USING THERMOPLASTICS**

PRIOR TO USE: Refer to instructions provided with each thermoplastic mask.


**UPDATE OVERLAY THERMOPLASTIC CAPABILITY**

1. To use Type-S™ thermoplastic, ensure threaded Silverman pins are installed in overlay and Posicast® locating buttons are removed from overlay.
2. To use Posicast® thermoplastic, ensure Posicast® locating buttons are installed in overlay and threaded Silverman pins are removed from overlay.

**THREADED SILVERMAN PINS**

1. To install or uninstall threaded Silverman pins, remove overlay from tabletop.
2. Locate threaded Silverman pins on underside of overlay and install or remove using flat head screwdriver.


3. Reattach overlay to tabletop.

**ATTACHING OVERLAY TO TABLETOP**

PRIOR TO USE: Before using this product, refer to operating instructions provided with MRI scanner.

- NOTE:
- If head coil support is not attached to overlay, see Attaching Head Coil Support section.
  - Tools Required: 3/16" Allen Wrench.


1. Remove table pads and fixed head coil from table.
2. Attach overlay supports to overlay.


3. Set overlay on table in desired location, ensure overlay supports rest in slots on side of table.

**WARNING**

- Grasp overlay at head end and foot end, using two people, to prevent pinched fingers.

4. If necessary, install shims to reduce lateral movement.

**WARNING**

- Avoid pinch points when overlay is in motion by ensuring patient does not grasp overlay sides.

**COIL SUPPORTS**

PRIOR TO USE: Refer to operating instructions provided with each RF Coil.


**ADJUSTING HEAD COIL SUPPORT**

1. Lift head coil support out of patient area by loosening adjustment knobs.
2. Small coil clip and large coil clip are interchangeable. If required, remove small coil clip and attach large coil clip.
  - a. Loosen top adjustment knob.
  - b. Remove small coil clip and bottom clamp support.
  - c. Attach bottom clamp support and large coil clip.
  - d. Tighten top adjustment knob.


3. Position patient on overlay.
4. Adjust head coil support with adjustment knobs and adjustment levers for desired coil placement.

NOTE:

- Ensure head coil support is secure.
- Head/Neck/Shoulder Type-S™ style thermoplastic mask may interfere with adjustment levers when Head Coil Supports are between indexing marks 22 and 32 and lower.


5. If necessary, anterior coil connector can be used to ensure coils connect when wrapped around patient.


**ATTACHING BODY COIL SUPPORT**

1. Position patient on overlay.
- NOTE: Ensure body coil support bridges are at highest vertical position.
2. Place body coil support bridges on overlay in desired position. Secure docking levers.
  3. Adjust height of body coil support bridges as needed. Secure bridge height with height adjustment knobs.


4. Adjust angular position of body coil support bridges, as needed. Secure bridge angular position with adjustment knobs.


**ATTACHING HEAD COIL SUPPORT TO OVERLAY**

1. Align head coil support piece with overlay rod.
2. Loosen adjustment lever on head coil support.
3. Carefully slide head coil support from head coil support piece onto overlay rod.

4. Discard head coil support piece accordingly.
5. Secure adjustment lever on head coil support.
6. Repeat steps 1 through 5 on opposite side.
7. Attach underneath head coil support.
  - a. Position support onto overlay rod.
  - b. Secure locking knob.

- c. Repeat steps a and b on opposite side.

NOTE: See *Attaching Overlay to Tabletop* section to proceed, including steps for initial setup.

**REPROCESSING**

**WARNING**  
 Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. If necessary, clean surface by removing visual contamination with common germicidal wipe, such as alcohol. If visual contamination cannot be removed repeat cleaning steps and if necessary, discard device.

**MAINTENANCE**

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2021 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. IPPS, LOK-BAR, POSICAST, PRODIGY AND TYPE-S ARE TRADEMARKS OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.