ProForm™ Overlay Patient Handles

INSTRUCTIONS FOR USE (IFU)

INTENDED USE

The device is indicated to aid in supporting, positioning, and/or immobilization of adult and pediatric patients undergoing radiation therapy of the head, brain, neck, and spine including radiosurgery and electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Do not use if device appears damaged.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Ensure device is secure prior to use.
- Ensure device/patient does not move for duration of set-up and treatment.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.


MRI SAFETY INFORMATION

MR
Conditional

- ProForm™ Overlay patient handle is MR Conditional. Device contains brass and must remain outside imaging region to prevent RF heating.

USING PROFORM™ OVERLAY PATIENT HANDLES

NOTE: Confirm ProForm™ Overlay patient handles are in correct position when changing from overlay to extension by consulting setup sheet.

1. Place ProForm™ Overlay patient handles onto overlay/extension at desired position.
2. Turn locking knob to secure handles in place.

REPROCESSING**⚠ WARNING**

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. If necessary, clean surface by removing visual contamination with common germicidal wipe, such as alcohol. If visual contamination cannot be removed repeat cleaning steps and if necessary, discard device.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. PROFORM IS A TRADEMARK OF MEDTEC LLC.
ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MEDTEC LLC