AccuForm™ Shim

REFERENCE GUIDE

CIVCO
Radiotherapy™

CE

INTENDED USE

The device is indicated to aid in supporting, positioning, and/or immobilization of adult and pediatric patients undergoing radiation therapy of the head, brain, neck, and spine including radiosurgery and electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

-  **WARNING**
- Do not use if device appears damaged.
 - Verification image should be taken to confirm patient position prior to treatment.
 - Ensure patient remains still for duration of setup and treatment.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

MRI SAFETY INFORMATION

- Device is MR safe.

USING ACCUFORM™ SHIM

1. Setup patient with AccuForm™ Cushion (*Refer to AccuForm™ Cushion instructions*).
2. Allow cushion to set (harden) around patient's head.

NOTE: Before placing shim, allow AccuForm™ Cushion to set (harden) for approximately 4 minutes.

3. Place AccuForm™ Shim between patient's head and AccuForm™ Cushion.

NOTE: Do not fold shim to increase thickness. If increased thickness is desired, stack 2 or more shims on top of each other.

4. Pull prepared thermoplastic mask on patient (*Refer to thermoplastic instructions*).
5. After allowing mask to cool, remove shim for all subsequent treatment.

REPROCESSING

-  **WARNING**
- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2022 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. CIVCO RADIOTHERAPY AND CIVCO RT ARE TRADEMARKS OF MEDTEC LLC. ACCUFORM IS A TRADEMARK OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.