# Respiratory Pump and Gauge

*Body Pro-Lok™ System*

---

REFERENCE GUIDE

**CIVCO**  
Radiotherapy™

CE

**INDICATIONS FOR USE**

The device is indicated to assist in the proper positioning of patients for radiation therapy including electron, photon and proton treatments, and simulation including CT and MR image acquisition.

**INTENDED USE**

The device is part of a system intended to immobilize, position and reposition patients undergoing stereotactic radiation therapy including SBRT.

**CAUTION**

Federal (United States) law restricts this device to sale by or on the order of a physician.

**⚠ WARNING**

- Do not use if device appears damaged.
- Ensure clearance is maintained between Body Pro-Lok™ system and all linac or imaging equipment before using system.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at [WWW.CIVCORT.COM](http://WWW.CIVCORT.COM).
- Verify patient position with completed setup sheet prior to treatment.
- Pressure indication is for reference only.
- Inspect device prior to use for signs of damage and general wear.

**NOTE:** If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

**MRI SAFETY INFORMATION**

MR  
Conditional

- Device is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:
  - Static magnetic field of 1.5 T and 3.0 T
  - Maximum spatial field gradient of 51,700 gauss/cm (517 T/m)

**BODY PRO-LOK™ RESPIRATORY BELT**

1. Attach respiratory gauge and pump to respiratory belt.
2. On gauge, use range A through H for respiratory belt inflation.
3. Inflate respiratory belt until desired respiratory restriction is achieved.
4. Prior to imaging or treatment, remove gauge and pump from respiratory belt.

**NOTE:** Note all indexed positions and acceptable pressure range on setup sheet.

**COMFORTCARE™ CUSHION**

1. Use range 1 through 8 on respiratory gauge for ComfortCare™ cushion inflation.
2. Inflate ComfortCare™ cushion to distribute pressure and improve patient comfort.
3. Prior to imaging or treatment, remove pump and gauge from cushion.

**NOTE:** Record all patient specific settings on setup sheet.

**REPROCESSING****⚠ WARNING**

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. Wipe all surfaces with common germicidal or antiseptic wipe, such as alcohol.

**MAINTENANCE**

**NOTE:** Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. CIVCO RADIOTHERAPY AND CIVCO RT ARE TRADEMARKS OF MEDTEC LLC. BODY PRO-LOK AND COMFORTCARE ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.