Evolve

INSTRUCTIONS FOR USE (IFU)

CE

INTENDED USE

The device is indicated to aid in supporting and positioning adult and adolescent patients undergoing radiation therapy of the breast and chest region including electron, photon, and proton treatments. Device is also used to position the patient during image acquisition to support treatment planning including in Computed Tomography (CT) and Magnetic Resonance (MR) imaging systems. The device is not intended for use with patients under 12 years of age.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Do not use if device appears damaged.
- Ensure tilting board is fixed properly before placing patient.
- Do not reposition device with patient on it.
- Do not exceed patient weight of 350 lbs (158.8 kg).
- Ensure device is secure prior to use.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Only use CQ Medical compatible accessories.
- Do not use excessive force on device. Injury or device damage may occur.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at www.CQmedical.com.
- Verify patient position with completed setup sheet prior to treatment.
- Do not allow patient to reposition themselves.
- Ensure angle adjustment is secure prior to handling.
- Ensure device is fully supported after adjusting angle.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

MRI SAFETY INFORMATION

- Device is MR safe.


MR Safe

ATTACHING EVOLVE™ SLANT BOARD

⚠ WARNING


- Lok-Bar™ is a positioning device only, do not use to secure patient.

1. Snap Lok-Bar™ at appropriate longitudinal position. Ensure Lok-Bar™ is secure.
2. Place Evolve™ onto Lok-Bar™ pins carefully to prevent damage (refer to Lok-Bar™ instruction). Ensure Evolve™ and Lok-Bar™ are secure.


USING EVOLVE™ SLANT BOARD

1. Twist clip manually to appropriate table width.


2. Lift Slant Board.
 3. Push down on lever to disengage and adjust to desired angle.
- NOTE:
- Ensure device is fully supported after adjusting angle.
 - Ensure Stabilization Bar rests flat on overlay/table.


USING EVOLVE™ BOTTOM BOARD

1. Secure Bottom Board to Slant Board using hex screws.


2. Lift and slide Bottom Stop to appropriate indexed location (0-14).

- NOTE:
- For patient loading/unloading, move Bottom Stop into inferior unregistered area. Ensure Bottom Stop is in a registered position at all other times.
 - Bottom Board can be folded.


3. To remove device from overlay/table, fold Bottom Board to Slant Board and lift. Ensure Bottom Board is secure during transport.

USING WING BOARDS WITH EVOLVE™

1. Place Wing Board on GlideTrack Lok-Bar™.


2. To adjust Wing Board, pull Wing Board Auto-Lok pin and slide to desired position.

- NOTE: If necessary, refer to specific Wing Board instruction.

FOR PRO WING BOARD (STEPS 3-4)

3. To adjust Pro Handles, pull handle Auto-Lok Pins and slide to desired position. Ensure pin is fully seated.

4. Place headrest at desired position using Wing Board attachment holes.

NOTE: Remove Wing Board from Slant Board before removing device from overlay.

REPROCESSING**⚠ WARNING**

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

1. If necessary, clean surface by removing visual contamination with 70% alcohol, Theracide (PDI Sani-Cloth Plus) wipe, PDI Super Sani Wipes, or PDI Bleach or PDI AF3 wipes. Thoroughly rinse with water after using wipes. If visual contamination cannot be removed, repeat cleaning steps and if necessary, discard device.

NOTE: Color may transfer to wipe during normal cleaning.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. LOK-BAR AND PRODIGY ARE TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS.
MEDTEC LLC