Bite Putty

For use with trUpoint ARCH™

INSTRUCTIONS FOR USE (IFU)

INTENDED USE

Device is intended to be used for immobilization, positioning and repositioning during Stereotactic Radiotherapy (SRT) and Stereotactic Radiosurgery (SRS) in all parts of the brain, head, and neck during external beam radiation therapy. Device is also used to immobilize and position the head during image acquisition to support treatment planning including in Computed Tomography (CT) and Magnetic Resonance (MR) imaging systems.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

⚠ WARNING

- Do not use if device appears damaged.
- Inspect single patient components for damage or contamination prior to use. If evident, dispose of as infectious waste and create new components.
- trUpoint ARCH™ is for use with CQ Medical specified accessories only.
- A confirmation image is recommended to ensure patient and target are in the desired location before treatment is delivered.

- NOTE:**
- If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.
 - For information regarding device accessories, please visit www.CQmedical.com.
 - Bite Cup (MTHF135), Bite Tray (MTHF033), Putty Dispenser (MTHF200)

MRI SAFETY INFORMATION

Device is MR safe.


MR Safe

- NOTE:** MR image quality may be compromised if the area of interest is relatively close to the bite cup and bite tray filled with putty, nasion cushion and arch. Optimize MR imaging parameters prior to treating patient. Image artifact extends approximately 19mm from bite cup and bite tray filled with putty, 16mm from nasion cushion and 12mm from the arch when scanned in nonclinical testing using a gradient echo sequence in a 3.0T Siemens Skyra MR system with built-in body coil.

FORMING BITE IMPRESSION TRAY**⚠ WARNING**

- Bite tray is intended for single patient throughout setup and treatment cycle.

1. Create solid dental impression using bite tray and putty (refer to putty instructions for additional warnings).
2. Ensure bite tray is completely full of putty.
3. Form bite impression in supine position.

**REPROCESSING****⚠ WARNING**

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.
- Inspect device for damage such as cracks or breakage. If damage is evident, discontinue use and contact CQ Medical for further instructions as needed.

BITE TRAY FILLED WITH PUTTY

1. Rinse with Listerine® or equivalent solution before and after use.

DISPOSAL**WARNING**

- *Dispose of single patient components as infectious waste.*

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COPYRIGHT © 2024 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. TRUPOINT ARCH IS A TRADEMARK OF MEDTEC LLC. POSICAST AND POSIFIX ARE REGISTERED TRADEMARKS OF MEDTEC LLC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
**MEDTEC LLC**