ThermFix™

INSTRUCTIONS FOR USE (IFU)

CE

Towards Gantry

Towards Gantry

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Kızağa Doğru

INTENDED USE

The device is indicated to aid in supporting and positioning adult and pediatric patients undergoing radiation therapy including electron, photon, and proton treatments. The device is also used during image acquisition to support treatment planning.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

 WARNING

- Do not use if device appears damaged.
- Ensure sliding clamps are retracted and covered or ThermFix™ inferior pelvic bar is removed prior to patient dismount.
- When positioning patient for first time, use setup sheet to record all adjustments. Setup sheet is available at www.CQmedical.com.
- Verify patient position with completed setup sheet prior to treatment.
- Verify all angles of treatment and attenuation characteristics prior to treating patients.
- Ensure device is secure prior to use.
- Ensure device/patient does not move for duration of set-up and treatment.
- Ensure gantry indicator is positioned towards gantry.
- Only use CQ Medical compatible accessories.
- Ensure accessories are securely attached in correct position prior to treatment and/or imaging.
- Do not allow patient to reposition themselves, such as arching back.

NOTE: If any serious incident occurs in relation to device, incident should be reported to manufacturer. If incident occurred within the European Union, also report to the competent authority of the Member State in which you are established.

MRI SAFETY INFORMATION


- Device is MR safe.


MR Safe

USING THERMIFIX™

1. Connect ThermFix™ superior bar marked "1" to ThermFix™ inferior bar marked "2" using ThermFix™ spacing tool.


USING THERMFIX™ SPACING TOOL


NOTE: ThermFix™ inferior pelvic bar can be attached to superior bar and docked to platform after patient is mounted on board.

2. Dock opposite side of ThermFix™ bars to UCT profile at desired location. Dock closer side of ThermFix™ bars by pulling tab and pushing down.

NOTE: Toward gantry arrow should face gantry.

3. Remove ThermFix™ spacing tool before pulling and/or docking masks.

4. Dock thermoplastic (see *thermoplastic instructions*).

- NOTE:
- Device is compatible with HipFix® thermoplastics.
 - To undock device; pull tab and lift.
 - ThermFix™ inferior pelvic bar can be removed before patient dismounts from board.
 - After initial setup, ThermFix™ spacing tool is not necessary. Place ThermFix™ superior and inferior bars at location recorded on setup sheet.

REPROCESSING

WARNING

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*

1. If necessary, clean surface by removing visual contamination with common germicidal wipe, such as alcohol. If visual contamination cannot be removed repeat cleaning steps and if necessary, discard device.

NOTE: Color may transfer to wipe during normal cleaning.

MAINTENANCE

NOTE: Inspect device prior to use for signs of damage and general wear.

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COPYRIGHT © 2023 ALL RIGHTS RESERVED. CQ MEDICAL IS A TRADEMARK OF MEDTEC LLC. THERMFX IS A TRADEMARK OF MEDTEC LLC.
ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
**MEDTEC LLC**